

Overview**Method Name**

Quantitative Enzymatic

NY State Available

Yes

Specimen**Specimen Type**

Varies

Specimen Required

Submit only 1 of the following specimens:

Serum

Draw blood in a plain, red-top tube(s). Serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Plasma

Draw blood in a green-top (sodium or lithium heparin) tube(s). Spin down and send 1 mL plasma refrigerated in a plastic vial.

Reject Due To

Hemolysis: Mild reject; Gross reject

Thawing: Warm reject; Cold OK

Lipemia: NA

Icterus: NA

Other: Body fluids

Specimen Minimum Volume

0.5 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	30 days	
	Frozen	30 days	
	Ambient	7 days	

Clinical & Interpretive**Reference Values**

Pancreatic amylase
6-35 months: 2-28 U/L
3-6 years: 8-34 U/L
7-17 years: 9-39 U/L
18 years and older: 12-52 U/L

Salivary amylase
18 months and older: 9-86 U/L

Total amylase
3-90 days: 0-30 U/L
3-6 months: 7-40 U/L
7-8 months: 7-57 U/L
9-11 months: 11-70 U/L
12-17 months: 11-79 U/L
18-35 months: 19-92 U/L
3-4 years: 26-106 U/L
5-12 years: 30-119 U/L
13 years and older: 30-110 U/L

Performance

PDF Report

No

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Test Classification

CPT Code Information

82150/x2