

## Overview

### Method Name

Quantitative Enzyme-Linked Immunosorbent Assay

### NY State Available

Yes

## Specimen

### Specimen Type

WB Heparin

### Specimen Required

Collect blood in a green top tube (sodium or lithium heparin). Submit 1 mL well-mixed blood in a plastic screw cap tube frozen.

**NOTE:** 1. Critical frozen. Separate samples must be submitted when multiple tests are ordered.

2. Unacceptable: non-frozen samples

### Reject Due To

Hemolysis: NA

Thawing: Warm reject; Cold reject

Lipemia: NA

Icterus: NA

Other: NA

### Specimen Minimum Volume

0.5 mL

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
WB Heparin	Frozen (preferred)	180 days	

**Clinical & Interpretive****Reference Values**

180 - 1800 nmol/L

**Performance****PDF Report**

No

**Performing Laboratory Location**

ARUP Laboratories

**Fees & Codes****Test Classification**

This test was developed and its performance characteristics determined by ARUP Laboratories. The U.S. Food and Drug Administration has not approved or cleared this test; however, FDA clearance or approval is not currently required for clinical use. The results are not intended to be used as the sole means for clinical diagnosis or patient management decisions.

**CPT Code Information**

83088

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
FHSTW	Histamine, Whole Blood	46436-2

Result ID	Reporting Name	LOINC®
Z2753	Histamine, Whole Blood	46436-2