

Overview**Method Name**

ImmunoCAP FEIA

NY State Available

Yes

Specimen**Specimen Type**

Serum

Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Reject Due To

Hemolysis	NA
Lipemia	Mild OK; Gross reject
Icterus	NA
Other	NA

Specimen Minimum Volume

0.5 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Clinical & Interpretive**Reference Values**

<0.35 kU/L

Interpretation

<u>Class</u>	<u>IgE (kU/L)</u>	<u>Comment</u>
0	<0.10	Negative
0/1	0.10 - 0.34	Equivocal
1	0.35 - 0.69	Low Positive
2	0.70 - 3.4	Moderate Positive
3	3.5 - 17.4	High Positive
4	17.5 - 49.9	Very High Positive
5	50.0 - 99.9	Very High Positive
6	> or = 100	Very High Positive

Performance**Method Description**

The test method is the Phadia ImmunoCAP allergen-specific IgE system.

PDF Report

No

Performing Laboratory Location

Eurofins Viracor

Fees & Codes**Test Classification**

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86003

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FDKYE	Dock Yellow (<i>Rumex crispus</i>) IgE	6738-9

Result ID	Reporting Name	LOINC®
Z3088	Dock Yellow (<i>Rumex crispus</i>) IgE	6738-9
Z3089	Class	15684-4