

## Overview

### Useful For

Monitoring therapy for kidney stones

Identifying increased urinary oxalate as a risk factor for stone formation

Diagnosis of primary or secondary hyperoxaluria

### Method Name

Only orderable as part of a profile. For more information see ROXUR / Oxalate, Random, Urine.

Enzymatic Using Oxalate Oxidase

### NY State Available

Yes

## Specimen

### Specimen Type

Urine

### Specimen Required

Only orderable as part of a profile. For more information see ROXUR / Oxalate, Random, Urine.

**Patient Preparation:** Avoid taking large doses (>2 g orally/24 hours) of vitamin C prior to specimen collection.

**Supplies:** Urine Tubes, 10 mL tube (T068)

**Container/Tube:** 10-mL plastic tube or a clean, plastic container with no metal cap

**Specimen Volume:** 7 mL

#### Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.
3. Specimen pH should be between 4.5 and 8 and will stay in this range if kept refrigerated. Specimens with pH above 8 may indicate bacterial contamination, and testing will be cancelled. Do not attempt to adjust pH as it will adversely affect results.

### Specimen Minimum Volume

1 mL

### Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	

**Clinical & Interpretive**
**Clinical Information**

Oxalate is an end product of glyoxalate and glycerate metabolism. Humans have no enzyme capable of degrading oxalate so it must be eliminated by the kidney.

In tubular fluid, oxalate can combine with calcium to form calcium oxalate stones. In addition, high concentrations of oxalate may be toxic for renal cells.

Increased urinary oxalate excretion results from inherited enzyme deficiencies (primary hyperoxaluria), gastrointestinal disorders associated with fat malabsorption (secondary hyperoxaluria), or increased oral intake of oxalate-rich foods or vitamin C.

Since increased urinary oxalate excretion promotes calcium oxalate stone formation, various strategies are employed to lower oxalate excretion.

**Reference Values**

Only orderable as part of a profile. For more information see ROXUR / Oxalate, Random, Urine.

No established reference values

**Interpretation**

An elevated urine oxalate (>0.46 mmol/day) may suggest disease states such as secondary hyperoxaluria (fat malabsorption), primary hyperoxaluria (alanine glyoxalate transferase enzyme deficiency, glyceric dehydrogenase deficiency), idiopathic hyperoxaluria, or excess dietary oxalate or vitamin C intake.

In stone-forming patients high urinary oxalate values, sometimes even in the upper limit of the normal range, are treated to reduce the risk of stone formation.

The urinary oxalate creatinine ratio varies widely in young children from <0.35 mmol/mL at birth, to <0.15 mmol/mL at 1 year, to <0.10 mmol/mL at 10 years, and <0.05 mmol/mL at 20 years of age (see table below).(1)

Oxalate/Creatinine (mg/mg)	
Age (year)	95th Percentile
0-0.5	<0.175
0.5-1	<0.139
1-2	<0.103
2-3	<0.08
3-5	<0.064
5-7	<0.056

7-17

&lt;0.048

**Cautions**

Ingestion of ascorbic acid (>2 g/day) may falsely elevate the measured urinary oxalate excretion.

**Clinical Reference**

1. Matos V, Van Melle G, Werner D, Bardy D, Guignard JP: Urinary oxalate and urate to creatinine ratios in a healthy pediatric population. *Am J Kidney Dis.* 1999;34:e1
2. Wilson DM, Liedtke RR: Modified enzyme-based colorimetric assay of urinary and plasma oxalate with improved sensitivity and no ascorbate interference: reference values and sample handling procedures. *Clin Chem.* 1991;37:1229-1235
3. Lieske JC, Wang X: Heritable traits that contribute to nephrolithiasis. *Urolithiasis.* 2019 Feb;47(1):5-10
4. Lieske JC, Turner ST, Edeh SN, Smith JA, Kardia SLR: Heritability of urinary traits that contribute to nephrolithiasis. *Clin J Am Soc Nephrol.* 2014 May;9(5):943-950. doi: 10.2215/CJN.08210813
5. Zhao F, Bergstralh EJ, Mehta RA, et al: Predictors of incident ESRD among patients with primary hyperoxaluria presenting prior to kidney failure. *Clin J Am Soc Nephrol.* 2016 Jan 7;11(1):119-126. doi: 10.2215/CJN.02810315

**Performance****Method Description**

The assay utilizes oxalate oxidase, which oxidizes oxalate to carbon dioxide and peroxide. In the presence of peroxidase, the peroxide oxidatively couples 3-methyl-2-benzothiazolinone and 3-dimethylaminobenzoic acid to form indamine dye, which is measured spectrophotometrically at 580 nm. (Kasidas GP, Rose GA: Continuous-flow assay for urinary oxalate using immobilized oxalate oxidase. *Ann Clin Biochem* 1985;22:412-419; package insert: Oxalate kit. Trinity Biotech; V 07/2016)

**PDF Report**

No

**Day(s) Performed**

Monday through Saturday

**Report Available**

3 days

**Specimen Retention Time**

7 days

**Performing Laboratory Location**

Rochester

**Fees & Codes**

**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been modified from the manufacturer's instructions. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.

**CPT Code Information**

83945

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
OXCO1	Oxalate, Random, U (mmol/L)	15086-2

Result ID	Test Result Name	Result LOINC® Value
OXCO1	Oxalate, Random, U (mmol/L)	15086-2