

## Overview

### Useful For

Identify women who are at high risk for developing preeclampsia.

### Method Name

Immunochemiluminometric Assay (ICMA)

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

Collection container/tube: 5 mL Red

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Acceptable: SST

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down within one hour and send 0.5 mL of serum frozen in a plastic vial.

### Reject Due To

Thawing\*\*                      Warm Reject; Cold Reject  
Other/Tissue/Swab          Specimens other than serum

### Specimen Minimum Volume

0.1 mL Note: This volume does not permit repeat analysis.

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Serum	Frozen (preferred)		
	Refrigerated		

## Clinical & Interpretive

### Reference Values

Age	Range (ng/mL)	
	Fasting	Random
Prepubertal	30-1000	10-500
Pubertal	20-200	20-100
Adults	10-150	0-40

## Performance

### PDF Report

No

### Performing Laboratory Location

Esoterix Endocrinology

## Fees & Codes

### Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

### CPT Code Information

83520

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FGBP1	IGF Binding Protein 1	12722-5

Result ID	Reporting Name	LOINC®
FGBP1	IGF Binding Protein 1	12722-5