

---

**Overview****Useful For**

An adjunct to urine D-lactate (preferred) for the diagnosis of D-lactate acidosis

**Special Instructions**

- [Biochemical Genetics Patient Information](#)

**Method Name**

Enzymatic

**NY State Available**

Yes

**Specimen****Specimen Type**

Plasma NaFI-KOx

**Ordering Guidance**

Urine is the preferred specimen for D-lactate determination, order DLAU / D-Lactate, Urine.

For determination of L-lactate (lactic acid), order LACS1 / Lactate, Plasma

**Specimen Required**

**Collection Container/Tube:** Sodium Fluoride/Potassium Oxalate Tube, 2 mL (T275)

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

**Collection Instructions:** Centrifuge, aliquot plasma into plastic vial, and freeze immediately.

**Reject Due To**

Gross hemolysis	OK
Gross lipemia	OK
Gross icterus	OK

---

**Specimen Minimum Volume**

0.55 mL

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Plasma NaFl-KOx	Frozen (preferred)	365 days	
	Ambient	7 days	
	Refrigerated	7 days	

**Clinical & Interpretive****Clinical Information**

D-lactate is produced by bacteria residing in the colon when carbohydrates are not completely absorbed in the small intestine. When large amounts of D-lactate are present, individuals can experience metabolic acidosis, altered mental status (from drowsiness to coma), and a variety of other neurologic symptoms, in particular dysarthria and ataxia. Although a temporal relationship has been described between elevations of plasma and urine D-lactate and the accompanying encephalopathy, the mechanism of neurologic manifestations has not been elucidated.

D-lactic acidosis is typically observed in patients with a malabsorptive disorder, such as short-bowel syndrome, or following a jejunioileal bypass. In addition, healthy children presenting with gastroenteritis may also develop the clinical presentation of D-lactic acidosis.

Routine lactic acid determinations in blood will not reveal abnormalities because most lactic acid assays measure only L-lactate. Accordingly, D-lactate analysis must be specifically requested (eg, this test). However, as D-lactate is readily excreted in urine, it is the preferred specimen for D-lactate determinations; see DLAU / D-Lactate, Urine.

**Reference Values**

0.0-0.25 mmol/L

**Interpretation**

Increased levels are consistent with D-lactic acidosis. However, because D-lactate is readily excreted, urine determinations are preferred.

**Cautions**

The test performed was D-lactate. This is a product of bacterial overgrowth in the gastrointestinal tract. It should not be confused with L-lactate, which accumulates in some metabolic acidosis.

---

**Clinical Reference**

1. Petersen C: [D-lactic acidosis. Nutr Clin Pract. 2005;20\(6\):634-645](#)

2. Kowlgi NG, Chhabra L: D-Lactic acidosis: An underrecognized complication of short bowel syndrome. *Gastroenterol Res Pract.* 2015;2015:476215. doi: /10.1155/2015/476215

**Performance****Method Description**

D-lactate is oxidized to pyruvate in the presence of D-lactate dehydrogenase and nicotinamide adenine dinucleotide (NAD<sup>+</sup>). The reaction proceeds because the pyruvate is continually removed as a pyruvate-hydrazone complex. The quantity of NADH produced is directly proportional to the amount of D-lactate oxidized and is measured spectrophotometrically at 340 nm.(Brandt RB, Siegel SA, Waters MG, Bloch MH: Spectrophotometric assay for D-(-)-lactate in plasma. *Anal Biochem.*1980;102(1):39-46; Cowan T, Pasquali M: Laboratory investigations of inborn errors of metabolism. In: K Sarafoglou, GF Hoffman, KS Roth, eds. *Pediatric Endocrinology and Inborn Errors of Metabolism.* 2nd ed. McGraw-Hill Education; 2017:1139-1158)

**PDF Report**

No

**Specimen Retention Time**

1 month

**Performing Laboratory Location**

Rochester

**Fees & Codes****Test Classification**

This test was developed, and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.

**CPT Code Information**

83605

**LOINC® Information**

---

Test ID	Test Order Name	Order LOINC Value
DLAC	D-Lactate, P	14045-9

Result ID	Reporting Name	LOINC®
8878	D-Lactate, P	14045-9