

## Overview

**Method Name**

Enzyme Immunoassay (EIA)

**NY State Available**

Yes

## Specimen

**Specimen Type**

Serum

**Specimen Required**

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 2 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s). (Serum gel tube is acceptable.) Spin down and send 2 mL serum refrigerate in a plastic vial

**Reject Due To**

Hem NA

olysis

Lipe NA

mia

Icter NA

us

Othe r Specimen that is too viscous to pipette. Tissue, biopsy, sputum, bronchial aspirate, FNA, bone marrow aspirate, stool or samples in transport media, fixative or Isolator tubes.

**Specimen Minimum Volume**

1.2 mL

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Serum	Refrigerated (preferred)	14 days	
	Ambient	14 days	
	Frozen		

## Clinical & Interpretive

### Reference Values

Reference interval: None Detected

Reportable Range: Positive Results reported in ng/mL from 0.20 ng/mL to 20.00 ng/mL.

Positive Results above 20.00 ng/mL are reported as "Above the Limit of Quantification".

### Cautions

Sodium hydroxide and sputolysin. Cross-reactivity occurs between blastomycosis and histoplasmosis, and in Paracoccidioidomycosis, penicilliosis, coccidioidomycosis, aspergillosis, and sporotrichosis.

## Performance

### PDF Report

No

### Performing Laboratory Location

MiraVista Diagnostics

## Fees & Codes

### Test Classification

This test was developed and its performance characteristics determined by MiraVista Diagnostics. It has not been cleared or approved by the FDA; however, FDA clearance or approval is not currently required for clinical use. The results are not intended to be used as the sole means for clinical diagnosis or patient management decisions.

### CPT Code Information

---

87385

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
FHIST	MVista Histoplasma Ag, S	51753-2

Result ID	Reporting Name	LOINC®
Z1711	Result:	51753-2
Z1035	Interpretation	Not Provided