

---

**Overview****Method Name**

Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

**NY State Available**

Yes

**Specimen****Specimen Type**

Varies

**Specimen Required****Submit only 1 of the following specimens:****Plasma****Specimen Type:** Plasma (Preferred)**Container/Tube:** Green-top (sodium heparin) tube(s).**Specimen volume:** 3 mL**Collection Instructions:** Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 3 mL sodium heparin plasma refrigerated in a plastic vial.**Serum****Specimen Type:** Serum**Container/Tube:** Red-top tube, serum gel is not acceptable.**Specimen volume:** 3 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 3 mL of serum refrigerated in a plastic vial.

**Reject Due To**

Hemolysis NA  
Lipemia NA  
Icterus NA  
Other NA

**Specimen Minimum Volume**

0.6 mL

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Frozen	180 days	
	Ambient		

**Clinical & Interpretive****Reference Values**

Units: ng/mL

Expected steady state concentrations of risperidone and 9-hydroxyrisperidone (combined total) in patients receiving recommended daily dosages: 10 - 120 ng/mL.

**Performance****PDF Report**

No

**Specimen Retention Time**

---

2 weeks

**Performing Laboratory Location**

Medtox Laboratories, Inc.

**Fees & Codes****Test Classification**

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

**CPT Code Information**

80342

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
FRISP	Risperidone and Metabolite	55157-2

Result ID	Reporting Name	LOINC®
Z1183	Risperidone (Risperdal)	9393-0
Z1056	9-Hydroxyrisperidone	9383-1
Z1057	Combined Total	9394-8