

**Overview****Method Name**

High Performance Liquid Chromatography with Fluorescence Detection (HPLC-FL)

**NY State Available**

Yes

**Specimen****Specimen Type**

Varies

**Specimen Required**

Submit only 1 of the following specimens:

**Plasma**

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

**Serum**

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

**Reject Due To**

Hemolysis NA

Lipemia NA

Icterus NA

Other NA

**Specimen Minimum Volume**

0.3 mL

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient		

**Clinical & Interpretive****Reference Values**

---

Reference Range: 500 - 4000 ng/mL

Serum Sotalol concentrations producing beta-blockade:

500 - 4000 ng/mL

Toxic range has not been established.

## Performance

### PDF Report

No

### Specimen Retention Time

2 weeks

### Performing Laboratory Location

Medtox Laboratories, Inc.

## Fees & Codes

### Test Classification

### CPT Code Information

80299

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FSOTA	Sotalol (Betapace)	12416-4

Result ID	Test Result Name	Result LOINC Value
Z1192	Sotalol	12416-4