

---

**Reporting Title:** Hereditary Prostate Cancer Panel**Performing Location** Rochester**Ordering Guidance:**

Customization of this panel and/or single gene analysis for any gene present on this panel is available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known mutations testing) is available for the genes on this panel. For more information see FMTT / Familial Mutation, Targeted Testing, Varies.

**Shipping Instructions:**

Specimen preferred to arrive within 96 hours of collection.

**Specimen Requirements:**

**Patient Preparation:** A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

**Specimen Type:** Whole blood

**Container/Tube:**

**Preferred:** Lavender top (EDTA) or yellow top (ACD)

**Acceptable:** Any anticoagulant

**Specimen Volume:** 3 mL

**Collection Instructions:**

1. Invert several times to mix blood.
2. Send specimen in original tube. **Do not** aliquot.

**Specimen Stability Information:** Ambient (preferred) 4 days/Refrigerated

**Forms:**

[1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available in Special Instructions:](#)

[-Informed Consent for Genetic Testing \(T576\)](#)

[-Informed Consent for Genetic Testing-Spanish \(T826\)](#)

2. [Molecular Genetics: Inherited Cancer Syndromes Patient Information](#)(T519) in Special Instructions

3. [Targeted Genes and Methodology Details for Hereditary Prostate Cancer Panel](#) in Special Instructions

4. If not ordering electronically, complete, print, and send a [Oncology Test Request](#) (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies (preferred)	0 hours	

**Supplemental Report:**

Supplemental

**CPT Code Information:**

81408

81162

81403

81292

81295

81298

81307

81319

81351

81479

**Reference Values:**

An interpretive report will be provided.