

Reporting Title: Transferrin, S  
Performing Location: Rochester

Specimen Requirements:  
Supplies: Sarstedt Aliquot Tube, 5 mL (T914)  
Container/Tube:  
Preferred: Serum gel  
Acceptable: Red top  
Submission Container/Tube: Plastic vial  
Specimen Volume: 1 mL

Collection Instructions:  
1. Within 2 hours of collection, centrifuge the specimen.  
2. For serum gel tubes, aliquot serum into a plastic vial prior to shipment.  
3. For red-top tubes, aliquot the serum into a plastic vial immediately after centrifuging.

Forms:  
If not ordering electronically, complete, print, and send a [Benign Hematology Test Request Form](#) (T755) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	7 days	
	Refrigerated (preferred)	7 days	
	Frozen	180 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
TRSF	Transferrin, S	Numeric	mg/dL	3034-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:  
No

CPT Code Information:  
84466

Reference Values:  
200-360 mg/dL