

**Reporting Title:** Ravulizumab, S  
**Performing Location:** Mayo Clinic Laboratories - Rochester Superior Drive

**Ordering Guidance:**  
To measure only serum concentration of ravulizumab, order RAVU / Ravulizumab, Serum.

To measure the impact of ravulizumab on complement activity and its effect on complement blockage, order RAVMP / Ravulizumab Monitoring Panel, Serum, which measures the alternative pathway function.

**Specimen Requirements:**  
**Patient Preparation:** Consider discontinuing natalizumab at least 4 weeks prior to specimen collection. Patient should consult the healthcare provider who prescribed this medication to determine if discontinuation is an option. If not, ok to proceed with testing while taking natalizumab.  
**Collection Container/Tube:**  
**Preferred:** Red top  
**Acceptable:** Serum gel  
**Submission Container/Tube:** Plastic vial  
**Specimen Volume:** 1 mL

- Collection Instructions:**
1. Draw blood immediately before next scheduled dose.
  2. Immediately after specimen collection, place the tube on wet ice.
  3. After specimen has clotted on wet ice, centrifuge at 4 degrees C and aliquot serum into a plastic vial.
  4. Freeze specimen within 30 minutes of centrifugation. Specimen must be placed on dry ice if not frozen immediately.

**Forms:**  
[If not ordering electronically, complete, print, and send a Coagulation Test Request](#) (T753) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	28 days	
	Ambient	28 days	
	Refrigerated	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
609420	Ravulizumab, S	Numeric	mcg/mL	97184-6

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
80299

**Reference Values:**

Lower limit of quantitation=5.0 mcg/mL

>175 mcg/mL-Therapeutic concentration for paroxysmal nocturnal hemoglobinuria and atypical hemolytic uremic syndrome