

Reporting Title: Urea, Random, U**Performing Location:** Mayo Clinic Laboratories - Rochester Main Campus**Ordering Guidance:**

A timed 24-hour urine collection is the preferred specimen for measuring and interpreting this urinary analyte. See URCR / Uric Acid, 24 Hour, Urine.

Random collections normalized to urinary creatinine may be of some clinical use in patients who cannot collect a 24-hour specimen, typically small children.

Specimen Requirements:**Supplies:** Sarstedt 5 mL Aliquot Tube (T914)**Container/Tube:** Plastic tube**Specimen Volume:** 4 mL**Collection Instructions:**

1. Collect a random urine specimen.
2. No preservative.

Forms:

If not ordering electronically, complete, print, and send a [Renal Diagnostics Test Request](#) (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Ambient	7 days	
	Frozen	30 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
URCON	Urea, Random, U	Numeric	mg/dL	3092-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84540

Reference Values:

No established reference values

Random urine urea may be interpreted in conjunction with serum urea, using both values to calculate fractional excretion of urea.

The calculation for fractional excretion (FE) of urea is

$$FE(U) = \frac{[U(urine) \times Creat(serum)]}{[U(serum) \times Creat(urine)]} \times 100$$