

**Reporting Title:** Adalimumab QN with Reflex to Ab, S**Performing Location:** Mayo Clinic Laboratories - Rochester Superior Drive**Ordering Guidance:**

If both quantitation and antibody testing are needed, regardless of the quantitation results, order ADALP / Adalimumab Quantitative with Antibody, Serum

**Specimen Requirements:**

**Patient Preparation:** For 12 hours before specimen collection, patient **should not** take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7)

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.5 mL Serum

**Collection Instructions:**

1. Draw blood immediately before the next dose of drug administration (trough specimen).
2. Centrifuge and aliquot serum into a plastic vial.

**Forms:**

If not ordering electronically, complete, print, and send one of the following with the specimen:

-[Gastroenterology and Hepatology Test Request](#) (T728)

-[Therapeutics Test Request](#) (T831)

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
ADALX	Adalimumab QN with Reflex to Ab, S	Numeric	mcg/mL	86894-3

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

80145

83520 (if appropriate)

**Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
ADLAB	Adalimumab Ab, S	1	83520	No	No

**Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
ADLAB	ADLAB	Adalimumab Ab, S	Numeric	AU/mL	90779-0

**Reference Values:**

ADALIMUMAB QUANTITATIVE:

Limit of quantitation is 0.8 mcg/mL. Optimal therapeutic ranges are disease specific.

ADALIMUMAB ANTIBODY:

&lt;14.0 AU/mL