

Reporting Title: Haloperidol, S
Performing Location: Mayo Clinic Laboratories - Rochester Superior Drive

Specimen Requirements:
Supplies: Sarstedt Aliquot Tube, 5 mL (T914)
Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions:

- 1. Draw blood immediately before next scheduled dose.
- 2. Within 2 hours of collection, centrifuge and aliquot serum into a plastic vial.

Forms:
If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
80339	Haloperidol, S	Numeric	ng/mL	3669-9
169	Reduced Haloperidol	Numeric	ng/mL	38364-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
80173

Reference Values:
Haloperidol:
5-17 ng/mL

Reduced Haloperidol:
10-80 ng/mL