

Reporting Title: Serotonin, S
Performing Location: Rochester

Additional Testing Requirements:
First-line testing for the diagnosis of carcinoid tumors with symptoms suggestive of carcinoid syndrome consists of urinary serotonin (SERU / Serotonin, 24 Hour, Urine), urinary 5-HIAA (HIAA / 5-Hydroxyindoleacetic Acid, 24 Hour, Urine), and serum chromogranin A (CGAK / Chromogranin A, Serum).

Specimen Requirements:
Patient Preparation:
1. Patients should not eat avocados, bananas, butternuts, cantaloupe, dates, eggplant, grapefruit, hickory nuts, honeydew melon, kiwifruit, melon, nuts, pecans, pineapple, plantains, plums, tomatoes, or walnuts, which are high in serotonin for 48 hours before and during collection.
2. Patient should discontinue medications that may elevate urine serotonin concentration including lithium, monoamine oxidase inhibitors, methyl dopa, morphine, and reserpine. Patient should discontinue use of selective serotonin reuptake inhibitors (SSRI eg, PROZAC), which can lead to depletion of platelet serotonin levels and result in false-negative serotonin results.
Supplies: Sarstedt Aliquot Tube, 5 mL (T914)
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 2.5 mL
Collection Instructions: Centrifuge as soon as blood has clotted and aliquot serum into plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	4 days	
	Refrigerated (preferred)	21 days	
	Frozen	90 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
84395	Serotonin, S	Numeric	ng/mL	27057-9

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
84260

Reference Values:

< or =230 ng/mL