

Test Definition: TP

Protein, Total, Serum

Reporting Title: Protein, Total, S

Performing Location: Mayo Clinic Laboratories - Rochester Main Campus

Necessary Information: Patient's age and sex are required.

Specimen Requirements:

Collection Container/Tube:
Preferred: Serum gel
Acceptable: Red top
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL
Collection Instructions:
1. Sorum gol tubos should be contrifuged

1. Serum gel tubes should be centrifuged within 2 hours of collection.

2. Red-top tubes should be centrifuged, and the serum aliquoted into a plastic vial within 2 hours of collection.

Forms:

If not ordering electronically, complete, print, and send a <u>Renal Diagnostics Test Request</u> (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	180 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
ТР	Protein, Total, S	Numeric	g/dL	2885-2

LOINC[®] and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information: 84155

Reference Values:

> or =1 year: 6.3-7.9 g/dL

Reference values have not been established for patients who are <12 months of age.