

Reporting Title: Chromium, S  
Performing Location: Rochester

**Ordering Guidance:**  
The US Food and Drug Administration recommended test for monitoring chromium in patients with metal-on-metal implants is CRWB / Chromium, Blood.

**Specimen Requirements:**  
**Patient Preparation:** High concentrations of gadolinium and iodine are known to interfere with most metal tests. If either gadolinium- or iodine-containing contrast media has been administered, a specimen should not be collected for 96 hours.  
**Supplies:** Metal Free Specimen Vial (T173)  
**Collection Container/Tube:** Plain, royal blue-top Vacutainer plastic trace element blood collection tube  
**Submission Container/Tube:** 7-mL Mayo metal-free, screw-capped, polypropylene vial  
**Specimen Volume:** 0.5 mL  
**Collection Instructions:**

1. Allow the specimen to clot for 30 minutes; then centrifuge the specimen to separate serum from the cellular fraction.
2. Remove the stopper. Carefully pour specimen into a Mayo metal-free, polypropylene vial, avoiding transfer of the cellular components of blood. **Do not** insert a pipet into the serum to accomplish transfer, and **do not** ream the specimen with a wooden stick to assist with serum transfer.
3. See [Metals Analysis Specimen Collection and Transport](#) for complete instructions.

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	28 days	METAL FREE
	Refrigerated (preferred)	28 days	METAL FREE
	Frozen	28 days	METAL FREE

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
8638	Chromium, S	Numeric	ng/mL	5622-6

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
82495

**Reference Values:**  
<0.3 ng/mL

When collected by a phlebotomist experienced in ultra-clean collection technique and handled according to the

---

instructions in [Metals Analysis Specimen Collection and Transport](#), we have observed the concentration of chromium in serum to be below 0.3 ng/mL. However, the majority of specimens submitted for analysis from unexposed individuals contain 0.3 ng/mL to 0.9 ng/mL of chromium. Commercial evacuated blood collection tubes not designed for trace-metal specimen collection yield serum containing 2.0 ng/mL to 5.0 ng/mL chromium derived from the collection tube.