



The following applies to Hereditary Cancer Custom Gene Panel. Testing is performed to evaluate for the presence of variants in coding regions and extending to +/- 10 base pairs of adjacent intronic sequence on either side of the coding exons of the genes analyzed. In addition, the analysis will cover select noncoding variants. Next-generation sequencing, multiplex ligation-dependent probe amplification (MLPA), and/or a polymerase chain reaction-based quantitative method is performed to test for the presence of copy number variants (CNV) in the genes analyzed. Polymerase chain reaction and gel electrophoresis is performed to test for the presence of the 10 megabase (Mb) inversion of coding exons 1–7 of the *MSH2* gene. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

This list is current from December 2024 to the present. This document is intended to highlight additional evaluations for variants of high clinical interest as well as technical limitations. However, this document does not comprehensively reflect all genomic regions covered by this test. For questions regarding transcripts, genes or regions covered, contact the laboratory at 800-533-1710.

Genomic Build: GRCh37 (hg19) unless otherwise specified

Gene	Reference Transcript	Additional Evaluations	Technical Limitations
<i>AIP</i>	NM_003977.4	-	-
<i>ALK</i>	NM_004304.5	-	-
<i>APC</i>	NM_000038.6	Promoter 1A: c.-172 to c.-19 (variants between c.-565 to c.-173 may be detected) Promoter 1B: c.-30632 to c.-30046	-
<i>ATM</i>	NM_000051.3	-	-
<i>AXIN2</i>	NM_004655.4	-	-
<i>BAP1</i>	NM_004656.4	-	-
<i>BARD1</i>	NM_000465.4	-	-
<i>BLM</i>	NM_000057.4	-	-
<i>BMPR1A</i>	NM_004329.2	-	-
<i>BRCA1</i>	NM_007294.4	+/- 20 base pairs of adjacent intronic sequence on either side of the coding exons	-
<i>BRCA2</i>	NM_000059.3	+/- 20 base pairs of adjacent intronic sequence on either side of the coding exons	-
<i>BRIP1</i>	NM_032043.3	-	-
<i>BUB1B</i>	NM_001211.5	-	-
<i>CDC73</i>	NM_024529.4	-	-
<i>CDH1</i>	NM_004360.5	-	-
<i>CDK4</i>	NM_000075.4	-	-
<i>CDKN1B</i>	NM_004064.4	-	-
<i>CDKN2A</i>	NM_000077.4	-	-
<i>CHEK2</i>	NM_007194.4	-	-
<i>CTNNA1</i>	NM_001903.5	-	-
<i>DICER1</i>	NM_177438.2	-	-
<i>DIS3L2</i>	NM_152383.4	-	Sequence variants and CNV in exon 19 may not be detected or reported
<i>EGFR</i>	NM_005228.5	-	-
<i>ELP1</i>	NM_003640.5	-	-

Targeted Genes and Methodology Details for Hereditary Cancer Custom Gene Panel (continued)

Gene	Reference Transcript	Additional Evaluations	Technical Limitations
<i>EPCAM</i>	NM_002354.3	-	Analysis for sequence variants will not be performed
<i>EXT1</i>	NM_000127.2	-	-
<i>EXT2</i>	NM_207122.1	-	-
<i>FANCA</i>	NM_000135.4	-	-
<i>FH</i>	NM_000143.3	-	-
<i>FLCN</i>	NM_144997.7	CNV in noncoding exons 1–3	-
<i>GPC3</i>	NM_004484.4	-	-
<i>GREM1</i>	NM_013372.7	-	Analyzed for CNV duplication of upstream enhancer region only; analyses for sequence variants and additional CNV will not be performed
<i>HOXB13</i>	NM_006361.5	-	-
<i>KIT</i>	NM_000222.2	-	-
<i>LZTR1</i>	NM_006767.4	-	-
<i>MAX</i>	NM_002382.5	-	-
<i>MC1R</i>	NM_002386.3	-	-
<i>MEN1</i>	NM_130799.2	-	-
<i>MET</i>	NM_001127500.3	-	-
<i>MITF</i>	NM_000248.3	-	Analyzed only for the presence of c.952G>A p.E318K (rs149617956); analysis for other sequence variants and CNV will not be performed
<i>MLH1</i>	NM_000249.3	c.-27C>A (rs587779001) c.454-13A>G (rs267607749) c.1990-16_1990-2del (rs267607881)	-
<i>MLH3</i>	NM_001040108.1	-	CNV in exon 5 will not be detected or reported
<i>MSH2</i>	NM_000251.3	10 Mb inversion of exons 1-7	-
<i>MSH3</i>	NM_002439.5	-	-
<i>MSH6</i>	NM_000179.2	-	-
<i>MUTYH</i>	NM_001128425.1	c.504+19_504+31del (rs781222233)	-
<i>NBN</i>	NM_002485.4	-	CNV in exon 16 will not be detected or reported
<i>NF1</i>	NM_000267.3	-	Analysis for sequence variants beyond coding exons +/- 10 base pairs of adjacent intronic sequence will not be performed
<i>NF2</i>	NM_000268.3	-	-
<i>NTHL1</i>	NM_002528.7	-	-
<i>PALB2</i>	NM_024675.4	-	-
<i>PDGFRA</i>	NM_006206.6	-	CNV in exons 17 and 19 will not be detected or reported

Targeted Genes and Methodology Details for Hereditary Cancer Custom Gene Panel (continued)

Gene	Reference Transcript	Additional Evaluations	Technical Limitations
<i>PHOX2B</i>	NM_003924.4	-	-
<i>PMS2</i>	NM_000535.7	-	-
<i>POLD1</i>	NM_002691.4	-	-
<i>POLE</i>	NM_006231.4	-	-
<i>POT1</i>	NM_015450.3	-	CNV in exon 5 will not be detected or reported
<i>PRKAR1A</i>	NM_002734.4	-	-
<i>PTCH1</i>	NM_000264.5	-	-
<i>PTEN</i>	NM_000314.8	Promoter: c.-1302 to c.-589	-
<i>RAD51B</i>	NM_133509.4	-	-
<i>RAD51C</i>	NM_058216.3	-	-
<i>RAD51D</i>	NM_002878.3	-	-
<i>RB1</i>	NM_000321.2	c.-240 to c.-1	-
<i>RECQL4</i>	NM_004260.3	-	-
<i>REST</i>	NM_005612.5	-	-
<i>RET</i>	NM_020975.6	-	-
<i>RNF43</i>	NM_017763.5	-	-
<i>SDHA</i>	NM_004168.4	-	-
<i>SDHAF2</i>	NM_017841.2	-	-
<i>SDHB</i>	NM_003000.3	-	-
<i>SDHC</i>	NM_003001.3	-	-
<i>SDHD</i>	NM_003002.4	-	-
<i>SMAD4</i>	NM_005359.6	-	-
<i>SMARCA4</i>	NM_001128849.2	-	CNV in exon 27 will not be detected or reported
<i>SMARCB1</i>	NM_003073.5	-	-
<i>SMARCE1</i>	NM_003079.5	-	CNV in exon 2 will not be detected or reported
<i>STK11</i>	NM_000455.5	-	-
<i>SUFU</i>	NM_016169.3	-	-
<i>TMEM127</i>	NM_017849.3	-	-
<i>TP53</i>	NM_000546.5	-	-
<i>TRIP13</i>	NM_004237.4	-	-
<i>TSC1</i>	NM_000368.5	-	-
<i>TSC2</i>	NM_000548.5	-	-
<i>VHL</i>	NM_000551.3	-	-
<i>WRN</i>	NM_000553.6	-	-
<i>WT1</i>	NM_024426.6	-	-

Targeted Genes and Methodology Details for Hereditary Cancer Custom Gene Panel (continued)

Available Hereditary Cancer Panels

Test ID	Test Name	Genes
BAP1Z	BAP1-Tumor Predisposition Syndrome, BAP1 Full Gene Analysis	<i>BAP1</i>
BHDZ	Birt-Hogg-Dube Syndrome, FLCN Full Gene Analysis	<i>FLCN</i>
BRGYP	Hereditary Breast/Gynecologic Cancer Panel	<i>ATM, BARD1, BRCA1, BRCA2, BRIP1, CDH1, CHEK2, EPCAM, MLH1, MSH2, MSH6, NBN, NF1, PALB2, PMS2, PTEN, RAD51C, RAD51D, STK11, TP53</i>
CDHZ	Hereditary Diffuse Gastric Cancer Syndrome, CDH1 Full Gene Analysis	<i>CDH1</i>
COMCP	Hereditary Common Cancer Panel	<i>APC, ATM, AXIN2, BARD1, BMPR1A, BRCA1, BRCA2, BRIP1, CDH1, CDK4, CDKN2A, CHEK2, DICER1, EPCAM, GREM1, HOXB13, MEN1, MLH1, MSH2, MSH3, MSH6, MUTYH, NBN, NF1, NTHL1, PALB2, PMS2, POLD1, POLE, PTEN, RAD51C, RAD51D, RET, SMAD4, STK11, TP53</i>
CRCGP	Hereditary Gastrointestinal Cancer Panel	<i>APC, ATM, AXIN2, BMPR1A, CDH1, CHEK2, CTNNA1, EPCAM, GREM1, KIT, MLH1, MLH3, MSH2, MSH3, MSH6, MUTYH, NTHL1, PDGFRA, PMS2, POLD1, POLE, PTEN, RNF43, SMAD4, STK11, TP53</i>
ENDCP	Hereditary Endocrine Cancer Panel	<i>AIP, APC, CDC73, CDKN1B, DICER1, FH, MAX, MEN1, NF1, PHOX2B, PRKAR1A, PTEN, RET, SDHA, SDHAF2, SDHB, SDHC, SDHD, TMEM127, TP53, TSC1, TSC2, VHL, WRN</i>
HBOCZ	BRCA1/BRCA2 Genes, Full Gene Analysis	<i>BRCA1, BRCA2</i>
HPGLP	Hereditary Paraganglioma/Pheochromocytoma Panel	<i>FH, MAX, NF1, RET, SDHA, SDHAF2, SDHB, SDHC, SDHD, TMEM127, VHL</i>
LRCCZ	Hereditary Leiomyomatosis and Renal Cell Cancer Syndrome, FH Full Gene Analysis	<i>FH</i>
LYNCP	Lynch Syndrome Panel	<i>MLH1, MSH2, MSH6, PMS2, EPCAM</i>
NF1Z	Neurofibromatosis Type 1, NF1 Full Gene Analysis	<i>NF1</i>
PANCP	Hereditary Pancreatic Cancer Panel	<i>ATM, BRCA1, BRCA2, CDKN2A, EPCAM, MLH1, MSH2, MSH6, PALB2, PMS2, STK11, TP53</i>
PRS8P	Hereditary Prostate Cancer Panel	<i>ATM, BRCA1, BRCA2, BRIP1, CHEK2, EPCAM, FANCA, HOXB13, MLH1, MSH2, MSH6, NBN, PALB2, PMS2, RAD51B, RAD51C, RAD51D, TP53</i>
PTNZ	PTEN Hamartoma Tumor Syndrome, PTEN Full Gene Analysis	<i>PTEN</i>
RENCP	Hereditary Renal Cancer Panel	<i>BAP1, DICER1, FH, FLCN, MET, MITF (c.952G>A p.E318K variant only), PTEN, SDHA, SDHAF2, SDHB, SDHC, SDHD, SMARCA4, SMARCB1, TMEM127, TP53, TSC1, TSC2, VHL</i>
RETZZ	Multiple Endocrine Neoplasia Type 2 Syndrome, RET Full Gene Analysis	<i>RET</i>
STK1Z	Peutz-Jeghers Syndrome, STK11 Full Gene Analysis	<i>STK11</i>
THYRP	Hereditary Thyroid Cancer Panel	<i>APC, DICER1, PRKAR1A, PTEN, RET, TP53, WRN</i>
VHLZZ	Von Hippel Lindau Syndrome, VHL Full Gene Analysis	<i>VHL</i>
WILMP	Hereditary Wilms Tumor Panel	<i>BLM, BUB1B, CDC73, DIS3L2, GPC3, REST, TP53, TRIP13, WT1</i>

Targeted Genes and Methodology Details

for Hereditary Cancer Custom Gene Panel (continued)

Available Hereditary Cancer Panels

Test ID	Test Name	Genes
XCP	Hereditary Expanded Cancer Panel	AIP, ALK, APC, ATM, AXIN2, BAP1, BARD1, BLM, BMPR1A, BRCA1, BRCA2, BRIP1, BUB1B, CDC73, CDH1, CDK4, CDKN1B, CDKN2A, CHEK2, CTTNA1, DICER1, DIS3L2, EGFR, ELP1, EPCAM, EXT1, EXT2, FANCA, FH, FLCN, GPC3, GREM1, HOXB13, KIT, LZTR1, MAX, MC1R, MEN1, MET, MITF (c.952G>A p.E318K variant only), MLH1, MLH3, MSH2, MSH3, MSH6, MUTYH, NBN, NF1, NF2, NTHL1, PALB2, PDGFRA, PHOX2B, PMS2, POLD1, POLE, POT1, PRKAR1A, PTCH1, PTEN, RAD51B, RAD51C, RAD51D, RB1, RECQL4, REST, RET, RNF43, SDHA, SDHAF2, SDHB, SDHC, SDHD, SMAD4, SMARCA4, SMARCEB1, SMARCE1, STK11, SUFU, TMEM127, TP53, TRIP13, TSC1, TSC2, VHL, WRN, WT1

Effective Date	Version	Synopsis of Test Change
June 2023	V2	Updated format (additional columns: Additional Evaluations and Technical Limitations)
December 2024	V3	FLCN: Added CNV detection of noncoding exons 1–3