
Overview

Useful For

Establishing or confirming the clinical diagnosis of hereditary hemochromatosis (HH) in adults

HFE genetic testing is NOT recommended for population screening

Testing of individuals with increased transferrin-iron saturation in serum and serum ferritin

With appropriate genetic counseling, predictive testing of individuals who have a family history of HH

Genetics Test Information

Detects the 2 common disease-causing mutations: C282Y and H63D. The S65C mutation is reported only when it is observed as part of the C282Y/S65C genotype.

Highlights

Molecular testing can be done to establish or confirm the diagnosis of hereditary hemochromatosis in individuals with clinical symptoms.

This test is not recommended for population screening.

This assay will not detect all of the mutations that cause hereditary hemochromatosis.

The S65C mutation is reported only when observed as part of the C282Y/S65C genotype.

Testing Algorithm

See [Hereditary Hemochromatosis Algorithm](#) in Special Instructions.

Special Instructions

- [Molecular Genetics: Congenital Inherited Diseases Patient Information](#)
- [Informed Consent for Genetic Testing](#)
- [Hereditary Hemochromatosis Algorithm](#)
- [Informed Consent for Genetic Testing \(Spanish\)](#)

Method Name

Polymerase Chain Reaction (PCR)-Based Assay Utilizing Agena Mass Array Platform

NY State Available

Yes

Specimen

Specimen Type

Varies

Shipping Instructions

Specimen preferred to arrive within 96 hours of draw.

Specimen Required

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant

Specimen Volume: 2.5 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send specimen in original tube.

Forms

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available in Special Instructions:

-[Informed Consent for Genetic Testing](#) (T576)

-[Informed Consent for Genetic Testing-Spanish](#) (T826)

2. [Molecular Genetics: Congenital Inherited Diseases Patient Information](#) (T521) in Special Instructions

3. If not ordering electronically, complete, print, and send a [Benign Hematology Test Request Form](#) (T755) with the specimen.

Specimen Minimum Volume

0.5 mL

Reject Due To

All specimens will be evaluated by Mayo Clinic Laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Frozen		
	Refrigerated		

For more information about hereditary hemochromatosis testing, see [Hereditary Hemochromatosis Algorithm](#) in Special Instructions.

Cautions

This assay will not detect all of the mutations that cause hereditary hemochromatosis. Therefore, the absence of a detectable mutation does not rule out the possibility that an individual is a carrier of or affected with this disease.

Test results should be interpreted in the context of clinical findings, family history, and other laboratory data. Errors in our interpretation of results may occur if information given is inaccurate or incomplete.

Rare polymorphisms exist that could lead to false-negative or false-positive results. If results obtained do not match the clinical findings, additional testing should be considered.

In rare cases, DNA alterations of undetermined significance may be identified.

Because of concerns of the overall penetrance of *HFE* mutations, *HFE* genetic testing is not recommended for population screening.

Clinical Reference

1. Mura C, Raguenes O, Ferec C: *HFE* Mutations analysis in 711 hemochromatosis probands: evidence for S65C implication in mild form of hemochromatosis. *Blood* 1999;93(8):2502-2505
2. Beutler E, Felitti VJ, Koziol J, et al: Penetrance of 845G->A (C282Y) *HFE* hereditary haemochromatosis mutation in the USA. *Lancet* 2002;359(9302):211-218
3. Walsh A, Dixon JL, Ramm GA, et al: The clinical relevance of compound heterozygosity for the C282Y and H63D substitutions in hemochromatosis. *Clin Gastroenterol Hepatol* 2006;4(11):1403-1410
4. Whitlock EP, Garlitz BA, Harris EL, et al: Screening for hereditary hemochromatosis: a systematic review for the U.S. Preventive Services Task Force. *Ann Intern Med* 2006;145(3):209-223

Performance

Method Description

[A PCR-based assay utilizing Agena Mass Array platform is used to test for the presence of C282Y, H63D, and S65C in the *HFE* gene. Because the S65C mutation has a minimal effect on iron metabolism, it is only reported when it is found with the C282Y mutation \(ie, if the patient has the C282Y/S65C genotype\). \(Unpublished Mayo method\)](#)

PDF Report

No

Day(s) and Time(s) Test Performed

Monday through Friday; 2 p.m.

Analytic Time

6 days

Maximum Laboratory Time

7 days

Specimen Retention Time

Whole Blood: 2 weeks (if available) Extracted DNA: 3 months

Performing Laboratory Location

Rochester

Fees and Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

81256-HFE (*hemochromatosis*) (eg, hereditary hemochromatosis) gene analysis, common variants (eg, C282Y, H63D)

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
HFE	Hemochromatosis HFE Gene Analysis, B	34519-9

Result ID	Test Result Name	Result LOINC Value
52899	Result Summary	50397-9
52900	Result	21694-5
52901	Interpretation	69047-9
52902	Specimen	31208-2
52903	Source	31208-2
52904	Method	49549-9
52905	Released By	18771-6