

Overview

Useful For

Assessing for *MYB* gene rearrangements in patients with primary salivary gland carcinoma to aid in confirming or excluding the diagnosis of primary salivary gland adenoid cystic carcinomas

Reflex Tests

Test ID	Reporting Name	Available Separately	Always Performed
_PBCT	Probe, +2	No, (Bill Only)	No
_PADD	Probe, +1	No, (Bill Only)	No
_PB02	Probe, +2	No, (Bill Only)	No
_PB03	Probe, +3	No, (Bill Only)	No
_IL25	Interphases,	No, (Bill Only)	No
_I099	Interphases, 25-99	No, (Bill Only)	No
_I300	Interphases, >=100	No, (Bill Only)	No

Testing Algorithm

This test does not include a pathology consultation. If a pathology consult is requested, PATHC / Pathology Consultation should be ordered and the appropriate FISH test will be ordered and performed at an additional charge.

This test includes a charge for application of the first probe set (2 FISH probes) and professional interpretation of results.

Â

Additional charges will be incurred for all reflex probes performed. Analysis charges will be incurred based on the number of cells analyzed per probe set. If no cells are available for analysis, no analysis charges will be incurred.

Method Name

Fluorescence In Situ Hybridization (FISH)

NY State Available

Yes

Specimen

Specimen Type

Tissue

Shipping Instructions

Advise Express Mail or equivalent if not on courier service.

Necessary Information

A reason for referral and pathology report are required in order for testing to be performed. Send information

with specimen. Acceptable pathology reports include working drafts, preliminary pathology or surgical pathology reports.

Specimen Required

Submit only 1 of the following specimens:

Specimen Type: Tissue

Preferred: Tissue block

Collection Instructions: Submit a formalin-fixed, paraffin-embedded (FFPE) tumor tissue block. Blocks prepared with alternative fixation methods may be acceptable; provide fixation method used.

Acceptable: Slides

Collection Instructions: Four consecutive, unstained, 5 micron-thick sections placed on positively charged slides, and 1 hematoxylin and eosin-stained slide.

Forms

If not ordering electronically, complete, print, and send an [Oncology Test Request](#) (T729) with the specimen.

Specimen Minimum Volume

Two consecutive unstained 5 micron-thick sections placed on positive-charged slides, and 1 hematoxylin and eosin-stained slide.

Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Tissue	Ambient (preferred)		
	Refrigerated		

Clinical and Interpretive

Clinical Information

Salivary adenoid cystic carcinomas (ACC), although uncommon, are frequent among salivary gland malignancies. ACC is typically an aggressive tumor with a poor prognosis. Histologically, ACC show significant morphologic overlap with other salivary gland tumors, but have a much different clinical course. Because ACC requires a management distinct from histologically similar lesions, it is important to make an accurate diagnosis. Translocations between *MYB* (6q23.3) and *NFIB* (9p24) have been identified in a large proportion of primary salivary gland ACC. These alterations have not been identified in other salivary gland tumors. Therefore, separation of *MYB*, in the proper clinical and histologic context, is diagnostic for ACC and can be confirmed by FISH with *MYB* break-apart probes.

Reference Values

An interpretive report will be provided.

Interpretation

A positive result is detected when the percent of cells with an abnormality exceeds the normal cutoff for the probe set.

A positive result suggests rearrangement of the *MYB* locus. The presence of a *MYB* rearrangement in conjunction with the proper clinical and histologic features is diagnostic of adenoid cystic carcinomas (ACC). A confirmed diagnosis of ACC results in specific clinical management that may be distinct from the management of other salivary gland neoplasms.

A negative result suggests no rearrangement of the *MYB* gene region at 6q23.3. The absence of a *MYB* rearrangement does not exclude the diagnosis of ACC, as a subset of ACCs do not show an *MYB* rearrangement.

Cautions

This test is not approved by the U.S. Food and Drug Administration, and it is best used as an adjunct to existing clinical and pathologic information.

This test is only intended to be used in the diagnosis of salivary gland tumors. This test has not been validated on adenoid cystic carcinoma (ACC) arising from nonsalivary gland locations. The results of this test are intended to be interpreted in association with the pathologic and clinical findings. The absence of *MYB* rearrangement does not exclude the diagnosis of ACC.

Fixatives other than formalin (eg, Prefer, Bouin) may not be successful for FISH assays, however nonformalin-fixed samples will not be rejected.

Paraffin-embedded tissues that have been decalcified are generally unsuccessful for FISH analysis. The pathologist reviewing the hematoxylin and eosin-stained slide may find it necessary to cancel testing.

Supportive Data

The probe was independently validated in a blinded study on 48 paraffin-embedded adenoid cystic carcinoma tissue specimens and 25 noncancerous control specimens. The normal controls were used to generate a normal cutoff for this assay.

Clinical Reference

1. Mitani Y, Rao PH, Futreal PA, et al: Novel chromosomal rearrangements and break points at the t(6;9) in salivary adenoid cystic carcinoma: association with *MYB-NFIB* chimeric fusion, MYB expression, and clinical outcome. *Clin Cancer Res* 2011;17(22):7003-7014
2. West RB, Kong C, Clarke N, et al: *MYB* expression and translocation in adenoid cystic carcinomas and other salivary gland tumors with clinicopathologic correlation. *Am J Surg Pathol* 2011;35(1):92-99
3. Bell D, Roberts D, Karpowicz M, et al: Clinical significance of Myb protein and downstream target genes in salivary adenoid cystic carcinoma. *Cancer Biol Ther* 2011;12:569-573

Performance

Method Description

This test is performed using a laboratory-developed *MYB* dual-color break-apart strategy probe (BAP). Paraffin-embedded tissues are cut at 5 microns and mounted on positively charged glass slides. The selection of tissue and

the identification of target areas on the hematoxylin and eosin (H and E)-stained slide are performed by a pathologist. Using the H and E-stained slide as a reference, target areas are etched with a diamond-tipped etcher on the back of the unstained slide to be assayed. The probe is hybridized to the appropriate target areas and 2 technologists each analyze 50 interphase nuclei (100 total) with the results expressed as the percent abnormal nuclei. (Unpublished Mayo method)

PDF Report

No

Day(s) and Time(s) Test Performed

Specimens are processed Monday through Sunday.

Results reported Monday through Friday, 8 a.m. to 5 p.m.

Analytic Time

7 days

Maximum Laboratory Time

10 days

Specimen Retention Time

Slides and H&E used for analysis are retained by the laboratory in accordance to CAP and NYS requirements. Client provided paraffin blocks and extra unstained slides (if provided) will be returned after testing is complete.

Performing Laboratory Location

Rochester

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed using an analyte specific reagent. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

88271x2, 88291 Æçâ,-â€œ DNA probe, each (first probe set), Interpretation and report

88271x2 Æçâ,-â€œ DNA probe, each; each additional probe set (if appropriate)

88271x1 Æçâ,-â€œ DNA probe, each; coverage for sets containing 3 probes (if appropriate)

88271x2 Æçâ,-â€œ DNA probe, each; coverage for sets containing 4 probes (if appropriate)

88271x3 Æçâ,-â€œ DNA probe, each; coverage for sets containing 5 probes (if appropriate)

88274 w/modifier 52 Æçâ,-â€œ Interphase in situ hybridization, <25 cells, each probe set (if appropriate)

88274 Interphase in situ hybridization, 25 to 99 cells, each probe set (if appropriate)

88275 Interphase in situ hybridization, 100 to 300 cells, each probe set (if appropriate)

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
SGTF	MYB (6q23), FISH, Ts	In Process

Result ID	Test Result Name	Result LOINC Value
54615	Result Summary	50397-9
54618	Interpretation	69965-2
54617	Result	62356-1
CG898	Reason for Referral	42349-1
54619	Specimen	31208-2
54620	Source	31208-2
54621	Tissue ID	80398-1
54622	Method	49549-9
55127	Additional Information	48767-8
53817	Disclaimer	62364-5
54623	Released By	18771-6