

Overview**Useful For**

Monitoring therapy

Determining compliance

Assessing toxicity

Method Name

Homogeneous Microparticle Agglutination Immunoassay

NY State Available

Yes

Specimen**Specimen Type**

Serum Red

Specimen Required**Container/Tube:** Red top**Specimen Volume:** 0.5 mL**Collection Instructions:** Tubes should be centrifuged and aliquoted within 2 hours of collection.**Forms**

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

[-Neurology Specialty Testing Client Test Request \(T732\)](#)[-Therapeutics Test Request \(T831\)](#)**Specimen Minimum Volume**

0.25 mL

Reject Due To

Gross hemolysis	Reject
-----------------	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	28 days	
	Ambient	48 hours	

Clinical and Interpretive

Clinical Information

Carbamazepine (Tegretol) is used in the control of partial seizures with both temporal lobe and psychomotor symptoms, and for generalized tonic-clonic seizures. It is also used for analgesia in trigeminal neuralgia.

Carbamazepine exhibits a volume of distribution of 1.4 L/kg with an elimination half-life of 15 hours. Protein binding averages 70%.

Carbamazepine-10,11-epoxide (CBZ10-11) is an active metabolite that represents the predominant form of the drug in children. The volume of distribution of CBZ10-11 is 1.1 L/kg, and the half-life is 5-8 hours.

Aplastic anemia and agranulocytosis are rare side effects of treatment with carbamazepine; baseline hematologic data should be documented before treatment is initiated.

Toxicity associated with carbamazepine overdose occurs when the blood level is 15.0 mcg/mL or higher and is typified by irregular breathing, muscle irritability, and hyperreflexia; followed by hyporeflexia, tachycardia, hypotension, and impaired consciousness with coma in severe toxicity. The higher the blood level, the more severe the symptoms.

Reference Values

Therapeutic: 4.0-12.0 mcg/mL

Critical value: > or =15.0 mcg/mL

Interpretation

Dosage adjustments are usually guided by monitoring blood levels. Most patients respond well when the serum concentration is in the range of 4.0 to 12.0 mcg/mL. Toxicity often occurs when levels are greater than or equal to 15.0 mcg/mL.

Cautions

Carbamazepine-10,11-epoxide (CBZ10-11) is not reported. A separate orderable test for this active metabolite is available (CARBG / Carbamazepine-10,11-Epoxide, Serum). Optimal response occurs when the CBZ10-11 serum level is in the range of 0.4 to 4.0 mcg/mL.

Clinical Reference

1. Cereghino JJ, Meter JC, Brock JT, et al: Preliminary observations of serum carbamazepine concentration in epileptic patients. *Neurology* 1973;23:357-366
2. Patsalos PN, Berry DJ, Bourgeois BF, et al: Antiepileptic drugs-best practice guidelines for therapeutic drug monitoring: A position paper by the subcommission on therapeutic drug monitoring, ILAE Commission on Therapeutic Strategies. *Epilepsia* 2008;49(7):1239-1276
3. Scheuer ML, Pedley TA: The evaluation and treatment of seizures. *N Engl JMed* 1990;322(21):1468-1474

Performance

Method Description

The ONLINE TDM Carbamazepine Gen.4 assay is a homogeneous microparticle agglutination immunoassay. It is a

2-reagent system used for the detection of carbamazepine in serum. Kinetic interaction of microparticles (KIMS) will be measured using automated analyzers. In this technology biotinylated drug hapten attached to streptavidin coated latex beads serves as the binding partner to anti-carbamazepine antibody. A competitive reaction to a limited amount of specific anti-carbamazepine antibody takes place between the latex bound hapten and free carbamazepine in the serum sample. A decrease in the apparent signal is proportional to the amount of drug present in the sample. (Package insert: Roche Carbamazepine reagent, Roche Diagnostic Corp, Indianapolis, IN 2016-04)

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

Same day/1 day

Specimen Retention Time

7 days

Performing Laboratory Location

Rochester

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test has been cleared, approved, or is exempt by the US Food and Drug Administration and is used per manufacturer's instructions. Performance characteristics were verified by Mayo Clinic in a manner consistent with CLIA requirements.

CPT Code Information

80156

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
CARTA	Carbamazepine, Tot, S	3432-2

Result ID	Test Result Name	Result LOINC Value
CARTA	Carbamazepine, Tot, S	3432-2