

Overview**Method Name**

Clot Dissolution

NY State Available

Yes

Specimen**Specimen Type**

Plasma Na Cit

Specimen Required

Collect blood in 3.2% sodium citrate light blue top tube. Centrifuge within 30 minutes after collection, freeze immediately. Send 2 mL platelet-poor plasma in plastic vial frozen.

Note:

1. Prohibit exercise prior to drawing sample.
2. To avoid release of plasminogen activator, do not massage vein vigorously, pump fist excessively or leave tourniquet in place for a prolonged period.

Specimen Minimum Volume

1 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	21 days	

Clinical and Interpretive**Reference Values**

>60 min

Performance**PDF Report**

No

Day(s) and Time(s) Test Performed

Wednesday

Analytic Time

3 days

Maximum Laboratory Time

5 - 18 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

CPT Code Information

85360

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FECLT	Euglobulin Clot Lysis Time	3244-1

Result ID	Test Result Name	Result LOINC Value
FECLT	Euglobulin Clot Lysis Time	3244-1