

**Overview**
**Method Name**

Qualitative Enzyme-Linked Immunosorbent Assay

**NY State Available**

Yes

**Specimen**
**Specimen Type**

Fecal

**Specimen Required**
**Preferred Specimen Type:** Unpreserved stool

**Supplies:** Sterile stool container

**Container/Tube:** Sterile stool container

**Specimen Volume:** 5 g

**Specimen Stability Information:** Refrigerated

**Collection Instructions:** 5 grams fresh, unpreserved stool or stool preserved in Cary-Blair transport media (Agar Swab is not acceptable), shipped refrigerate in a plastic leak-proof container.

**Specimen Minimum Volume**

1 gm

**Reject Due To**

Hemolysis:	NA
Thawing:	Cold OK; Warm reject
Lipemia:	NA
Icterus:	NA
Other:	Preservative other than Cary-Blair

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Fecal	Refrigerated (preferred)	14 days	
	Ambient	14 days	
	Frozen	14 days	

## Clinical and Interpretive

### Reference Values

Negative

### Interpretation

A positive result is indicative of the presence of lactoferrin, a marker for fecal leukocytes.

A negative result does not exclude the presence of intestinal inflammation.

## Performance

### PDF Report

No

### Day(s) and Time(s) Test Performed

Sunday through Saturday

### Analytic Time

1 - 2 days

### Maximum Laboratory Time

3 - 6 days

### Performing Laboratory Location

ARUP Laboratories

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### CPT Code Information

83630

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FLACF	Lactoferrin, Fecal by ELISA	40703-1

Result ID	Test Result Name	Result LOINC Value
FLACF	Lactoferrin, Fecal by ELISA	40703-1