

## Overview

### Useful For

Diagnosing congenital alpha-2 plasmin inhibitor deficiencies (rare)

Providing a more complete assessment of disseminated intravascular coagulation, intravascular coagulation and fibrinolysis, or hyperfibrinolysis (primary fibrinolysis), when measured in conjunction with fibrinogen, fibrin D-dimer, fibrin degradation products, soluble fibrin monomer complex, and plasminogen

Evaluating liver disease

Evaluating the effects of fibrinolytic or antifibrinolytic therapy

### Special Instructions

- [Coagulation Guidelines for Specimen Handling and Processing](#)

### Method Name

Chromogenic

### NY State Available

Yes

## Specimen

### Specimen Type

Plasma Na Cit

### Specimen Required

See [Coagulation Guidelines for Specimen Handling and Processing](#) in Special Instructions.

**Specimen Type:** Platelet-poor plasma

**Collection Container/Tube:** Light-blue top (citrate)

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

### Collection Instructions:

1. Centrifuge, remove plasma, and centrifuge plasma again.
2. Freeze specimen immediately at < or = -40 degrees C, if possible.

### Additional Information:

1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
2. Each coagulation assay requested should have its own vial.

**Forms**

If not ordering electronically, complete, print, and send a [Coagulation Test Request](#) (T753) with the specimen.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

|                 |        |
|-----------------|--------|
| Gross hemolysis | Reject |
| Gross lipemia   | Reject |
| Gross icterus   | Reject |

**Specimen Stability Information**

| Specimen Type | Temperature | Time    | Special Container |
|---------------|-------------|---------|-------------------|
| Plasma Na Cit | Frozen      | 14 days |                   |

**Clinical and Interpretive**
**Clinical Information**

Alpha-2 plasmin inhibitor (antiplasmin) is synthesized in the liver with a biological half-life of approximately 3 days. It inactivates plasmin, the primary fibrinolytic enzyme responsible for remodeling the fibrin thrombus, and binds fibrin, together with factor XIIIa, making the clot more difficult to lyse. Absence of alpha-2 plasmin inhibitor results in uncontrolled plasmin-mediated breakdown of the fibrin clot and is associated with increased risk of bleeding.

**Reference Values**

Adults: 80-140%

Normal, full-term, and premature infants may have mildly decreased levels (> or =50%) which reach adult levels < or = 90 days postnatal.\*

\*See Pediatric Hemostasis References section in [Coagulation Guidelines for Specimen Handling and Processing](#) in Special Instructions.

**Interpretation**

Patients with congenital homozygous deficiency (with levels of <10%) are clinically affected (bleeding). Heterozygotes having levels of 30% to 60% of mean normal activity are usually asymptomatic.

Lower than normal levels may be suggestive of consumption due to activation of plasminogen and its inhibition by alpha-2 plasmin inhibitor.

The clinical significance of high levels of alpha-2 plasmin inhibitor is unknown.

**Cautions**

Alpha-2 plasmin inhibitor results are potentially affected by:

-Heparin, unfractionated or low-molecular-weight >4 U/mL

-Alpha-2-macroglobulin >7 mg/mL; potentially leading to a falsely-increased result

-Hemoglobin >200 mg/dL

-Bilirubin >20 mg/dL

-Triglycerides >1000 mg/dL

### Clinical Reference

1. Lijnen HR, Collen D: Congenital and acquired deficiencies of components of the fibrinolytic system and their relation to bleeding or thrombosis. *Blood Coagul Fibrinolysis* 1989;3:67-77

2. Francis RB Jr: Clinical disorders of fibrinolysis: A critical review. *Blut* 1989;59:1-14

3. Aoki N: Hemostasis associated with abnormalities of fibrinolysis. *Blood Rev* 1989;3:11-17

### Performance

### Method Description

This assay is performed using the HemosIL Plasmin Inhibitor Kit on the Instrumentation Laboratory ACL TOP. Patient plasma, containing alpha-2 plasmin inhibitor, is mixed with reagent containing excess plasmin. Plasmin activity in the reagent is rapidly inhibited by alpha-2 plasmin inhibitor. Residual plasmin activity is then measured using an amidolytic activity assay, in which residual plasmin lyses a synthetic chromogenic substrate and subsequently releases paranitroaniline (detected at 405 nm) to a level that is inversely proportional to the amount of alpha-2 plasmin inhibitor in the sample. (Teger-Nilsson AC, Friberger P, Gyzander E: Determination of a new rapid plasmin inhibitor in human blood by means of a plasmin specific tripeptide substrate. *Scand J Clin Lab Invest* 1977;37:403-409)

### PDF Report

No

### Day(s) and Time(s) Test Performed

Monday through Friday; Varies

### Analytic Time

3 days

### Maximum Laboratory Time

7 days

### Specimen Retention Time

7 days

### Performing Laboratory Location

Rochester

### Fees and Codes

#### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.

- 
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
  - Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been modified from the manufacturer's instructions. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

85410

**LOINC® Information**

| Test ID | Test Order Name              | Order LOINC Value |
|---------|------------------------------|-------------------|
| A2PI    | Alpha-2 Plasmin Inhibitor, P | 27810-1           |

| Result ID | Test Result Name             | Result LOINC Value |
|-----------|------------------------------|--------------------|
| A2PI      | Alpha-2 Plasmin Inhibitor, P | 27810-1            |