

ZBTB46 IHC, Tech Only

### Overview

### **Useful For**

Identification of dendritic cells and classification of histiocytic disorders

#### Reflex Tests

Test ID	Reporting Name	Available Separately	Always Performed
IHTOA	IHC Additional, Tech Only	No	No
IHTOI	IHC Initial, Tech Only	No	No

### **Method Name**

Immunohistochemistry

### **NY State Available**

Yes

## **Specimen**

## **Specimen Type**

**TECHONLY** 

### **Advisory Information**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

## **Shipping Instructions**

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

## Specimen Required

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope

**Preferred:** 2 unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block

## **Digital Image Access**

- 1. Information on accessing digital images of IHC stains and the manual requisition form can be accessed through this website: <a href="https://www.mayocliniclabs.com/test-info/ihc/index.html">www.mayocliniclabs.com/test-info/ihc/index.html</a>
- 2. Clients ordering stains using a manual requisition form will not have access to digital images.



ZBTB46 IHC, Tech Only

3. Clients wishing to access digital images must place the order for IHC stains electronically. Information regarding digital imaging can be accessed through this website: <a href="https://www.mayocliniclabs.com/test-info/ihc/faq.html">www.mayocliniclabs.com/test-info/ihc/faq.html</a>

### **Forms**

If not ordering electronically, complete, print, and send an <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH)</u> <u>Stains Request</u> (T763) with the specimen.

## **Reject Due To**

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	
	Cytology smears
	Nonformalin fixed tissue
	Nonparaffin embedded tissue
	Noncharged slides
	ProbeOn slides
	Wet/frozen tissue

## **Specimen Stability Information**

Specimen Type	Temperature	Time
TECHONLY	Ambient (preferred)	
	Refrigerated	

## **Clinical and Interpretive**

### **Clinical Information**

Zinc finger and BTB domain-containing protein 46 (ZBTB46) is a transcription factor expressed by classical dendritic cells and committed dendritic cell precursors. Expression of ZBTB46 differentiates histiocytic disorders such as Langerhans cell histiocytosis, histiocytic sarcoma, and indeterminate cell histiocytosis from Rosai-Dorfman and Erdheim-Chester disease.

## Interpretation

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

The positive and negative controls are verified as showing appropriate immunoreactivity. If a control tissue is not included on the slide, a scanned image of the relevant quality control tissue is available upon request. Contact 855-516-8404.



ZBTB46 IHC, Tech Only

Interpretation of this test should be performed in the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### **Cautions**

No significant cautionary statements

### **Clinical Reference**

- 1. Reizis B: Classical dendritic cells as a unique immune cell lineage. J Exp Med 2012 209(6):1053-1056
- 2. Satpathy AT, Wumesh KC, Albring JC, et al: Zbtb46 expression distinguishes classical dendritic cells and their committed progenitors from other immune lineages. J Exp Med 2012 209(6):1135-1152
- 3. Meredith MM, Liu K, Darrasse-Jeze G, et al: Expression of the zinc finger transcription factor zDC (Zbtb46, Btbd4) defines the classical dendritic cell lineage. J Exp Med 2012 209(6):1156-1165
- 4. Satpathy AT, Brown RA, Gomulia E, et al: Expression of the transcription factor ZBTB46 distinguishes human histiocytic disorders of classical dendritic cell origin. Mod Path 2018 31:1479-1486

### **Performance**

### **Method Description**

Immunohistochemistry on sections of paraffin-embedded tissue.(Unpublished Mayo method)

## **PDF Report**

No

## Day(s) and Time(s) Test Performed

Monday through Friday

### **Analytic Time**

1 day

## **Maximum Laboratory Time**

3 days

### **Specimen Retention Time**

Until staining is complete

## **Performing Laboratory Location**

Rochester

### **Fees and Codes**

### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact Customer Service.

#### **Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with



ZBTB46 IHC, Tech Only

CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

## **CPT Code Information**

88342-TC, primary

88341-TC, if additional IHC

## **LOINC®** Information

Test ID	Test Order Name	Order LOINC Value
BTB46	ZBTB46 IHC, Tech Only	Order only;no result

Result ID	Test Result Name	Result LOINC Value
605255	ZBTB46 IHC, Tech Only	Bill only; no result