

## Overview

### Useful For

Measuring dabigatran concentration in plasma

This assay is **not useful for** measurement of other direct thrombin inhibitors eg, argatroban or bivalirudin.

### Special Instructions

- [Coagulation Guidelines for Specimen Handling and Processing](#)

### Method Name

Chromogenic Assay

### NY State Available

Yes

## Specimen

### Specimen Type

Plasma Na Cit

### Advisory Information

This test measures dabigatran only. For measurement of argatroban direct thrombin inhibitor, see ARGAT / Argatroban, Ecarin, Plasma.

This test is **not indicated for** monitoring low molecular weight heparin (LMWH), unfractionated heparin (UFH), or oral direct anti-Xa inhibitors (eg, apixaban, rivaroxaban, edoxaban). For monitoring oral direct anti-Xa inhibitors, see APIXA / Apixaban, Anti-Xa, Plasma; EDOXA / Edoxaban, Anti-Xa, Plasma, or RIVAR / Rivaroxaban, Anti-Xa, Plasma.

### Necessary Information

If priority specimen, mark request form, give reason, and request a call-back.

### Specimen Required

**Specimen Type:** Platelet-poor plasma

**Collection Container/Tube:** Light-blue top (3.2% sodium citrate)

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

### Collection Instructions:

1. Specimen should be drawn 1 to 3 hours (peak) after a dose or just prior (trough) to the next dose for dabigatran concentrations.

2. For complete instructions, see [Coagulation Guidelines for Specimen Handling and Processing](#) in Special Instructions.

3. Centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again.
4. Aliquot plasma into a plastic vial leaving 0.25 mL in the bottom of centrifuged vial.
5. Freeze plasma immediately (no longer than 4 hours after collection) at -20 degrees C or, ideally, < or =-40 degrees C.

**Additional Information:**

1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
2. Each coagulation assay requested should have its own vial.

**Forms**

If not ordering electronically, complete, print, and send a [Coagulation Test Request](#) (T753) with the specimen.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	42 days	

**Clinical and Interpretive**
**Clinical Information**

Dabigatran, an oral direct thrombin inhibitor (DTI) that directly inhibits factor IIa, is indicated for use in patients with nonvalvular atrial fibrillation and venous thromboembolism. Dabigatran is administered orally twice daily, is eliminated primarily through the renal system, and can inhibit both soluble and clot-bound thrombin. Dabigatran does not require routine therapeutic monitoring. However, in selected clinical situations, measurement of drug level would be useful, eg, renal insufficiency, assessment of compliance, periprocedural, suspected overdose, advanced age, and extremes of body weight.

**Observed dabigatran steady-state exposure concentrations**

Nonvalvular atrial fibrillation(1)	Median trough levels	Median peak levels
	(10-16 hours post dose)	(1-3 hours post dose)
110 mg twice daily	64.7 ng/mL (28.2-155)	126 ng/mL (52-275)
150 mg twice daily	91 ng/mL (39.8-215)	175 ng/mL (74.3-383)

<b>Venous thromboembolism(2)</b>	<b>Median trough levels</b> (10-16 hours post dose)	
150 mg twice daily	59.7 ng/mL (38.6-146)	

### Reference Values

<20 ng/mL

### Interpretation

Therapeutic reference ranges have not been established. See Clinical Information for peak and trough drug concentrations observed in clinical trials of dabigatran.

### Cautions

Dabigatran concentration may be affected by drug interactions as well as liver and renal disease.

Marked presence of lipemia or bilirubin in the sample could falsely decrease dabigatran levels.

### Supportive Data

The lower limit of detection of this assay is 20 ng/mL.

### Clinical Reference

1. Reilly PA, Lehr T, Haertter S, et al: The effect of dabigatran plasma concentrations and patient characteristics on the frequency of ischemic stroke and major bleeding in atrial fibrillation patients: the RE-LY Trial. *Am J Coll Cardiol.* 2014;63:321-328
2. Pradaxa (dabigatran etexilate). Public Assessment Record. European Medicines Agency; 2008 Update January 23, 2020. Accessed June 19, 2020. Available at [www.ema.europa.eu/en/medicines/human/EPAR/pradaxa](http://www.ema.europa.eu/en/medicines/human/EPAR/pradaxa)
3. Van Cott EM, Roberts AJ, Dager WE: Laboratory monitoring of parenteral direct thrombin inhibitors. *Semin Thromb Hemost.* 2017;43: 270-276
4. Gosselin RC, King JH, Janatpour KA, et al: Comparing direct thrombin inhibitors using aPTT, Ecarin clotting times, and thrombin inhibitor management testing. *Ann Pharmacother.* 2004;38:1383-1388
5. Gosselin RC, Adcock DM, Bates SM, et al: International Council for Standardization in Haematology (ICSH) recommendations for laboratory measurement of direct oral anticoagulants. *Thromb Haemost.* 2018 Mar;118(3):437-450
6. Lind SE, Boyle ME, Fisher S, Ishimoto J, Trujillo TC, Kiser TH: Comparison of the aPTT with alternative tests for monitoring direct trombin inhibitors in patient samples. *Am J Clin Pathol.* May 2014;141:665-674
7. Curvers J, van de Kerkhof D, Stroobants AK, van den Dool EJ, Scharnhorst V: Measuring direct thrombin inhibitors with routine and dedicated coagulation assays. *Am J Clin Pathol.* 2012;138: 551-558
8. Fenyvesi T, Jorg I, Harenberg J: Monitoring of anticoagulant effects of direct thrombin inhibitors. *Semin Thromb Hemost.* 2002;28(4):361-368

9. PRADAXA (dabigatran etexilate). Package insert: Boehringer Ingelheim Pharmaceuticals, Inc; 2019

## Performance

### Method Description

The dabigatran, Ecarin chromogenic assay is performed on the Instrumentation Laboratory ACL TOP 700 using the Diagnostica Stago ECA II kit. The STA ECA II kit is a chromogenic assay based on the cleavage of prothrombin by ecarin to meizothrombin which then enzymatically cleaves the paranitroaniline (pNA) from the chromogenic substrate resulting in a measureable colorometric change.(Package insert: STA -CA II. Diagnostica Stago S.A.S; Revision 09/2015)

### PDF Report

No

### Day(s) and Time(s) Test Performed

Monday through Friday; Varies

### Analytic Time

1 day

### Maximum Laboratory Time

3 days

### Specimen Retention Time

7 days

### Performing Laboratory Location

Rochester

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

### CPT Code Information

80299

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
DABIE	Dabigatran, Ecarin, P	In Process



---

Result ID	Test Result Name	Result LOINC Value
DABI1	Dabigatran, Ecarin, P	74220-5
DABI2	Interpretation	69049-5
DABI3	Cautions	62364-5