

## Overview

### Useful For

Detection of drug abuse involving phencyclidine (street names: angel dust, hog, or angel hair) in urine specimens handled through the chain-of-custody process

### Additional Tests

Test ID	Reporting Name	Available Separately	Always Performed
COCH	Chain of Custody Processing	No	Yes
ADLTX	Adulterants Survey, CoC, U	Yes	Yes

## Testing Algorithm

Testing for adulterants will be performed on all chain of custody urine samples per regulatory requirements.

### Method Name

Immunoassay/Gas Chromatography-Mass Spectrometry (GC-MS) Confirmation with Quantitation

### NY State Available

Yes

## Specimen

### Specimen Type

Urine

### Advisory Information

This test is for situations that require the chain-of-custody process. For testing **not** requiring chain of custody, order PCPU / Phencyclidine Confirmation, Random, Urine.

### Specimen Required

**Supplies:** Chain of Custody Kit (T282)

**Container/Tube:** Chain-of-Custody Kit containing the specimen containers, seals, and documentation required.

**Specimen Volume:** 10 mL

**Collection Instructions:** Collect specimen in the container provided, seal, and submit with the associated documentation to satisfy the legal requirements for chain-of-custody testing.

### Forms

1. [Chain-of-Custody Request](#) is included in the Chain-of-Custody Kit (T282).
2. If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

**Specimen Minimum Volume**

5 mL

**Reject Due To**

Gross hemolysis	OK
Gross icterus	OK

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	72 hours	

**Clinical and Interpretive**
**Clinical Information**

Phencyclidine (PCP) is a drug of abuse. This compound affects diverse neural pathways and interacts with cholinergic, adrenergic, GABA-secreting, serotonergic, opiate neuronal receptors, and gamma receptors. It has analgesic, anesthetic, and stimulatory effects, yielding bizarre behavior, ranging from depression through catatonia, euphoria, violent rage, and hallucinations. Most fatalities result from its hypertensive effect.

Diagnosis of PCP usage depends on drug screening. PCP is excreted in the urine.

Chain of custody is required whenever the results of testing could be used in a court of law. Its purpose is to protect the rights of the individual contributing the specimen by demonstrating that it was under the control of personnel involved with testing the specimen at all times; this control implies that the opportunity for specimen tampering would be limited. This includes a record of the disposition of a specimen to document the personnel who collected it, who handled it, and who performed the analysis. When a specimen is submitted in this manner, analysis will be performed in such a way that it will withstand regular court scrutiny.

**Reference Values**

Negative

Cutoff concentrations

IMMUNOASSAY SCREEN

&lt;25 ng/mL

PHENCYCLIDINE BY GC-MS

&lt;10 ng/mL

**Interpretation**

The presence of phencyclidine (PCP) in urine is a strong indicator that the patient has used PCP.

**Cautions**

Urine phencyclidine (PCP) may be undetectable at alkaline pH. Urine pH must, therefore, always be recorded.

**Clinical Reference**

1. Schuster DI, Arnold FJ, Murphy RB: Purification, pharmacological characterization and photoaffinity labeling of sigma receptors from rat and bovine brain. Brain Res 1995;670:14-28
2. Bayorh MA, Zokowska-Grojec A, Palkovits M, Kopin IJ: Effect of phencyclidine (PCP) on blood pressure and catecholamine levels in discrete brain nuclei. Brain Res 1984;321:315-318
3. Baselt RC: Disposition of Toxic Drugs and Chemicals in Man. 11th edition. Biomedical Publications, 2017

**Performance****Method Description**

The specimen is initially screened by immunoassay followed by gas chromatography/mass spectrometry.(Unpublished Mayo method)

**PDF Report**

No

**Day(s) and Time(s) Test Performed**

Wednesday

**Analytic Time**

2 days

**Maximum Laboratory Time**

8 days

**Specimen Retention Time**

14 days

**Performing Laboratory Location**

Rochester

**Fees and Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

83992

G0480 (if appropriate)

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
PCPX	Phencyclidine Confirmation, CoC, U	16254-5

Result ID	Test Result Name	Result LOINC Value
6672	Phencyclidine Immunoassay Screen	19659-2
36229	Phencyclidine-by GC/MS	16254-5
36230	Phencyclidine Interpretation	69050-3
36231	Chain of Custody	77202-0