

Overview**Method Name**

ImmunoCAP FEIA

NY State Available

Yes

Specimen**Specimen Type**

Serum

Specimen Required

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Clinical and Interpretive**Clinical Information**

This assay is used to detect allergen specific-IgE using the ImmunoCAP FEIA method. In vitro allergy testing is the

primary testing mode for allergy diagnosis.

Reference Values

Class IgE	(kU/L)	Comment
0	<0.10	Negative
0/1	0.10-0.34	Equivocal/Borderline
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	High Positive
4	17.50-49.99	Very High Positive
5	50.00-99.99	Very High Positive
6	>99.99	Very High Positive

Performance

Method Description

The ImmunoCAP FEIA method uses as the solid phase a flexible, hydrophobic cellulosic polymer to which allergen has been covalently linked. The advantage of this system is that it has a very high antigen binding capacity when compared to other systems and it has minimal non-specific binding with high total IgE. Viracor Eurofins provides an optional low range calibrator at 0.1 kU/L and a 0/1 class. This test has been cleared or approved for diagnostic use by the U.S. Food and Drug Administration.

PDF Report

No

Day(s) and Time(s) Test Performed

Monday through Friday

Analytic Time

1-2 days

Maximum Laboratory Time

4-7 days

Performing Laboratory Location

Viracor Eurofins

Fees and Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared

or approved by the U.S. Food and Drug Administration.

CPT Code Information

86008 x 4

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FHZCP	Hazelnut Component Panel	Not Provided

Result ID	Test Result Name	Result LOINC Value
Z5641	Hazelnut Component rCor a 1	69421-6
Z5642	Class	81995-3
Z5643	Hazelnut Component rCor a 8	58753-5
Z5644	Class	81993-8
Z5645	Hazelnut Component rCor a 9	65765-0
Z5646	Class	82002-7
Z5647	Hazelnut Component rCor a 14	81788-2
Z5648	Class	82569-5