

Overview

Useful For

Diagnosis of arsenic intoxication in 24-hour urine specimens

Testing Algorithm

Total arsenic will be performed first. If the total arsenic concentration is 15 mcg/L or greater, then fractionation will be performed and reported. If total arsenic is below 15 mcg/L, total arsenic will be reported as less than 15 mcg/L and fractionation will not be performed.

Special Instructions

- [Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens](#)
- [Trace Metals Analysis Specimen Collection and Transport](#)

Method Name

Ion Chromatography Extraction/Inductively Coupled Plasma-Mass Spectrometry (IC-ICP-MS)

NY State Available

Yes

Specimen

Specimen Type

Urine

Necessary Information

24 Hour volume is required.

Specimen Required

Patient Preparation:

1. Patient should not eat seafood for a 48-hour period prior to start of collection.
2. High concentrations of gadolinium and iodine are known to interfere with most metals tests. If either gadolinium- or iodine-containing contrast media has been administered, a specimen should not be collected for 96 hours.

Supplies: Urine Tubes, 10 mL (T068)

Collection Container/Tube: Clean, plastic urine container with no metal cap or glued insert

Submission Container/Tube: Plastic, 10-mL urine tube (T068) or clean, plastic aliquot container with no metal cap or glued insert

Specimen Volume: 10 mL

Collection Instructions:

1. Collect urine for 24 hours.

2. Refrigerate specimen within 4 hours of completion of 24-hour collection.
3. See [Trace Metals Analysis Specimen Collection and Transport](#) in Special Instructions for complete instructions.

Additional Information: See [Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens](#) in Special Instructions for multiple collections.

Urine Preservative Collection Options

Note: The addition of preservative or application of temperature controls **must occur within 4 hours of completion** of the collection.

Ambient	OK
Refrigerate	Preferred
Frozen	OK
50% Acetic Acid	OK
Boric Acid	No
Diazolidinyl Urea	No
6M Hydrochloric Acid	OK
6M Nitric Acid	OK
Sodium Carbonate	No
Thymol	No
Toluene	No

Specimen Minimum Volume

3 mL

Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	72 hours	

Clinical and Interpretive

Clinical Information

Arsenic (As) exists in a number of different forms; some are toxic, while others are not. The toxic forms are the

inorganic species of As(+3) (As-III), As(+5) (As-V), and their partially detoxified metabolites, monomethylarsine and dimethylarsine. As-III is more toxic than As-V and both are more toxic than mono- and dimethylarsine. The biologic half-life of inorganic arsenic is 4 to 6 hours, while the biologic half-life of the methylated metabolites is 20 to 30 hours. Target organs of As-III-induced effects are the heart, gastrointestinal tract, skin and other epithelial tissues, kidney, and nervous system.

Inorganic arsenic is carcinogenic to humans. Symptoms of chronic poisoning, called arseniasis, are mostly insidious and nonspecific. The gastrointestinal tract, skin, and central nervous system are usually involved. Nausea, epigastric pain, colic abdominal pain, diarrhea, and paresthesias of the hands and feet can occur.

Nontoxic, organic forms of arsenic are present in many foods. Arsenobetaine and arsenocholine are the 2 most common forms of organic arsenic found in food. The most common foods that contain significant concentrations of organic arsenic are shellfish and other predators in the seafood chain (cod, haddock, etc). Some meats such as chicken that have been fed on seafood remnants may also contain the organic forms of arsenic.

Following ingestion of arsenobetaine and arsenocholine, these compounds undergo rapid renal clearance to become concentrated in the urine. Organic arsenic is completely excreted within 1 to 2 days after ingestion and there are no residual toxic metabolites. The biologic half-life of organic arsenic is 4 to 6 hours.

For reporting purposes, the concentrations of the inorganic forms (As[+3] and As[+5]) along with the methylated forms (monomethylarsine and dimethylarsine) will be summed and reported together as "Inorganic" arsenic. This is consistent with how the biological exposure index (BEI) reference range is reported.

Reference Values

TOTAL ARSENIC

<18 mcg/24 hour

INORGANIC ARSENIC

<20 mcg/24 hour

Reference values apply to all ages.

*Biological exposure indices (BEI) for arsenic is 35 mcg/L based on the concentration of inorganic arsenic plus methylated metabolites.

Interpretation

The quantitative reference range for fractionated arsenic applies only to the inorganic forms. Concentrations greater than or equal to 20 mcg inorganic arsenic per 24 hours are considered toxic.

There is no limit to the normal range for the organic forms of arsenic, since they are not toxic and normally present after consumption of certain food types. For example, a typical finding in a urine specimen with total 24-hour excretion of arsenic of 350 mcg/24 hours would be that more than 95% is present as the organic species from a dietary source, and less than 5% present as the inorganic species. This would be interpreted as indicating the elevated total arsenic was due to ingestion of the nontoxic form of arsenic, usually found in food.

A normal value for blood arsenic does not exclude a finding of an elevated urine inorganic arsenic, due to the very short half-life of blood arsenic.

Cautions

Consumption of seafood before collection of a urine specimen for arsenic testing is likely to result in a report of an

elevated concentration of arsenic, which can be clinically misleading.

Clinical Reference

Caldwell KL, Jones RL, Verdon CP, et al: Levels of urinary total and speciated arsenic in the US population: National Health and Nutrition Examination Survey 2003-2004. J Expo Sci Environ Epidemiol 2009;19:59-68

Performance

Method Description

Inorganic (AsIII, AsV), organic (arsenobetaine, arsenocholine), and methylated (monomethylarsonic acid, dimethylarsinic acid) arsenic species in urine are separated using ion chromatography (IC) and quantitated by inductively coupled plasma-mass spectrometry (ICP-MS) in kinetic energy discrimination sensitive (KEDS) mode using an internal standard and a mobile phase matrix calibration. (Nixon DE, Moyer TP: Routine clinical determination of lead, arsenic, cadmium, and thallium in urine and whole blood by inductively coupled plasma mass spectrometry. Spectrochimica Acta B 1996;51:13-25; Sen I, ZOU W, Alvaran J, et al: Development and validation of a simple and robust method for arsenic speciation in human urine using HPLC/ICP-MS. J AOAC Int 2015;98[2]:517-523; Hata A, Endo Y, Nakajima Y, et al: HPLC-ICP-MS speciation analysis of arsenic in urine of Japanese subjects without occupational exposure. J Occup Health 2007;49[3]:217-223; Yoshimura Y, Endo Y, Shimoda Y, et al: Acute arsine poisoning confirmed by speciation analysis of arsenic compounds in the plasma and urine by HPLC-ICP-MS. J Occup Health 2011;53[1]:45-49)

PDF Report

No

Day(s) and Time(s) Test Performed

Monday through Friday; 8 a.m.

Analytic Time

2 days

Maximum Laboratory Time

4 days

Specimen Retention Time

14 days

Performing Laboratory Location

Rochester

Fees and Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

82175

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
ASFR	Arsenic Fractionation, 24 Hr, U	In Process

Result ID	Test Result Name	Result LOINC Value
3628	Inorganic Arsenic (Toxic)	53838-9
3629	Organic Arsenic (Non-Toxic)	53839-7
113133	Total Arsenic	5587-1
TM4	Collection Duration	13362-9
VL2	Urine Volume	3167-4