

Overview

Useful For

Detection of fucosidosis

This test is **not useful for** establishing carrier status for fucosidosis.

Genetics Test Information

Fucosidosis is an autosomal recessive lysosomal storage disorder caused by reduced or absent alpha-L-fucosidase enzyme activity.

Determining enzymatic activity is the next step of the diagnostic workup for an individual clinically suspicious for an oligosaccharidosis and with a positive screening result suggestive of fucosidosis.

Testing Algorithm

See [Lysosomal Storage Disorders Diagnostic Algorithm, Part 1](#) in Special Instructions.

Special Instructions

- [Informed Consent for Genetic Testing](#)
- [Biochemical Genetics Patient Information](#)
- [Informed Consent for Genetic Testing \(Spanish\)](#)
- [Lysosomal Storage Disorders Diagnostic Algorithm, Part 1](#)

Method Name

Fluorometric

NY State Available

Yes

Specimen

Specimen Type

Whole Blood ACD

Advisory Information

If clinically suspicious of an oligosaccharidosis, a screening test is available. Order OLIGU / Oligosaccharide Screen, Random, Urine.

Shipping Instructions

For optimal isolation of leukocytes, it is recommended the specimen arrive refrigerate within 144 hours of collection to be stabilized. Collect specimen Monday through Thursday only and not the day before a holiday. Specimen should be Collected and packaged as close to shipping time as possible.

Specimen Required

Container/Tube:

Preferred: Yellow top (ACD solution B)

Acceptable: Yellow top (ACD solution A)

Specimen Volume: 6 mL

Collection Instructions: Send specimen in original tube. Do not transfer blood to other containers.

Forms

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available in Special Instructions:

-[Informed Consent for Genetic Testing](#) (T576)

-[Informed Consent for Genetic Testing-Spanish](#) (T826)

2. [Biochemical Genetics Patient Information](#) (T602) in Special Instructions

3. If not ordering electronically, complete, print, and send an [Inborn Errors of Metabolism Test Request](#) (T798) with the specimen.

Specimen Minimum Volume

5 mL

Reject Due To

Gross hemolysis	Reject
-----------------	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Whole Blood ACD	Refrigerated (preferred)	6 days	YELLOW TOP/ACD
	Ambient	4 days	YELLOW TOP/ACD

Clinical and Interpretive

Clinical Information

Fucosidosis is an autosomal recessive lysosomal storage disorder caused by reduced or absent alpha-L-fucosidase enzyme activity. This enzyme is involved in degrading asparagine-linked, fucose-containing complex molecules (oligosaccharides and glycoasparagines) present in cells. Reduced or absent activity of this enzyme results in the abnormal accumulation of these molecules in the tissues and body fluids.

Severe and mild subgroups of fucosidosis, designated types I and II, have been described, although recent data suggests individual patients may represent a continuum within a wide spectrum of severity. The more severe type is characterized by infantile onset, rapid psychomotor regression, and severe neurologic deterioration. Additionally, dysostosis multiplex and elevated sweat sodium chloride are frequent findings. Death typically occurs within the first decade of life. Those with the milder phenotype express comparatively mild psychomotor and neurologic regression, radiologic signs of dysostosis multiplex and skin lesions (angiokeratoma corporis diffusum). Normal sweat salinity, the presence of the skin lesions, and survival into adulthood most readily distinguish milder from more severe phenotypes. Although the disorder is panethnic, the majority of reported patients with fucosidosis have been from Italy and southwestern United States. To date, about 100 cases have been reported worldwide.

An initial diagnostic workup includes a urine screening assay for several oligosaccharidosis (OLIGU / Oligosaccharide Screen, Random, Urine). If the screening assay is suggestive of fucosidosis, enzyme analysis of alpha-L-fucosidase can confirm the diagnosis.

Reference Values

> or =0.32 nmol/min/mg protein

Interpretation

Values below 0.32 nmol/min/mg protein are consistent with a diagnosis of fucosidosis.

Cautions

No significant cautionary statements

Clinical Reference

1. Cowan TM, Yu C: Laboratory investigations of inborn errors of metabolism. In Pediatric Endocrinology and Inborn Errors of Metabolism. Edited by K Sarafoglou, GF Hoffmann, KS Roth, New York, McGraw-Hill Medical Division, 2009, pp 867-868
2. Enns GM, Steiner RD, Cowan TM: Lysosomal disorders. In Pediatric Endocrinology and Inborn Errors of Metabolism. Edited by K Sarafoglou, GF Hoffmann, KS Roth, New York, McGraw-Hill Medical Division, 2009, pp 747-748
3. Thomas GH: Disorders of Glycoprotein Degradation: alpha-Mannosidosis, beta-Mannosidosis, Fucosidosis, and Sialidosis. In The Online Metabolic and Molecular Bases of Inherited Disease. Edited by D Valle, AL Beaudet, B Vogelstein, et al. New York, McGraw-Hill. Accessed March 18, 2019, Available at <https://ommbid.mhmedical.com/content.aspx?sectionid=225545029&bookid=2709&Resultclick=2>

Performance**Method Description**

Incubation of 4-methylumbelliferyl-alpha-L-fucopyranoside with cell homogenates results in cleavage of the substrate by alpha-L-fucosidase yielding 4-methylumbelliferone (4-MU) and fucose. Free 4-MU can be quantitated by measurement of the fluorescence.(Beratis NG, Turner BM, Labadie G, Hirschhorn K: a-L-fucosidase in cultured skin fibroblasts from normal subjects and fucosidosis patients. *Pediatr Res* 1977;11:862-866; Cowan T, Pasquali M: Laboratory Investigations of Inborn Errors of Metabolism. In Pediatric Endocrinology and Inborn Errors of Metabolism. Second edition. Edited by K Sarafoglou, GF Hoffman, KS Roth. 2017, pp 1139-1158)

PDF Report

No

Day(s) and Time(s) Test Performed

Varies

Analytic Time

30 days

Maximum Laboratory Time

45 days

Specimen Retention Time

WBC homogenate stored 1 month

Performing Laboratory Location

Rochester

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

82657

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FUCW	Alpha-Fucosidase, Leukocytes	24047-3

Result ID	Test Result Name	Result LOINC Value
8814	Alpha-Fucosidase, Leukocytes	24047-3
35635	Interpretation (FUCW)	59462-2
35634	Reviewed By	18771-6