

## Overview

### Useful For

Diagnosing deficiency of coagulation factor X, congenital or acquired

Evaluating hemostatic function in liver disease

Investigation of prolonged prothrombin time or activated partial thromboplastin time

### Special Instructions

- [Coagulation Guidelines for Specimen Handling and Processing](#)

### Method Name

Optical Clot-Based

### NY State Available

Yes

## Specimen

### Specimen Type

Plasma Na Cit

### Advisory Information

Coagulation testing is highly complex, often requiring the performance of multiple assays and correlation with clinical information. For that reason we suggest ordering Coagulation Consultations.

### Necessary Information

If priority specimen, mark request form, give reason, and request a call-back.

### Specimen Required

See [Coagulation Guidelines for Specimen Handling and Processing](#) in Special Instructions.

**Specimen Type:** Platelet-poor plasma

**Collection Container/Tube:** Light-blue top (citrate)

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

### Collection Instructions:

1. Within 4 hours of collection, centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again. Aliquot plasma into separate plastic vial leaving 0.25 mL in the bottom of centrifuged vial.

2. Freeze plasma immediately (no longer than 4 hours after collection) at -20 degrees C, or, ideally at < or = -40 degrees C.

**Additional Information:**

1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
2. Each coagulation assay requested should have its own vial.

**Forms**

[If not ordering electronically, complete, print, and send a Coagulation Test Request \(T753\)](#) with the specimen.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Gross hemolysis	Reject
Gross lipemia	Reject
Gross icterus	Reject

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

**Clinical and Interpretive**
**Clinical Information**

Factor X is a vitamin K-dependent serine protease that is synthesized in the liver. Its biological half-life is 24 to 48 hours. Factor X participates in both intrinsic and extrinsic pathways of coagulation (final common pathway) by serving as the enzyme (factor Xa) in the prothrombinase complex.

Congenital factor X deficiency is rare. Acquired deficiency associated with liver disease, warfarin therapy, vitamin K deficiency, systemic amyloidosis and inhibitors (rare). Deficiency may cause prolonged prothrombin time and activated partial thromboplastin time.

**Reference Values**

Adults: 70-150%

Normal, full-term newborn infants or healthy premature infants may have decreased levels (> or =15-20%) which may not reach adult levels for > or =180 days postnatal.\*

\*See Pediatric Hemostasis References section in [Coagulation Guidelines for Specimen Handling and Processing](#) in Special Instructions.

**Interpretation**

Acquired deficiency is more common than congenital deficiency.

Homozygotes: <25%

Heterozygotes: 25% to 50%

### Cautions

Liver disease, warfarin therapy, or vitamin K deficiency may lower factor X levels.

### Clinical Reference

1. Girolami A, Scandellari R, Scapin M, Vettore S: Congenital bleeding disorders of the vitamin K-dependent clotting factors. *Vitam Horm* 2008;78:281-374
2. Brenner B, Kuperman AA, Watzka M, Oldenburg J: Vitamin K-dependent coagulation factors deficiency. *Semin Thromb Hemost* 2009 Jun;35(4):439-446
3. Menegatti M, Peyvandi F: Factor X deficiency. *Semin Thromb Hemost* 2009 Jun;35(4):407-415
4. Girolami A, Ruzzon E, Tezza F, et al: Congenital FX deficiency combined with other clotting defects or with other abnormalities: a critical evaluation of the literature. *Haemophilia* 2008;14(2):323-328
5. Girolami A, Scarparo P, Scandellari R, Allemand E: Congenital factor X deficiencies with a defect only or predominantly in the extrinsic or in the intrinsic system: a critical evaluation. *Am J Hematol* 2008;83(8):668-671

### Performance

#### Method Description

[The factor X assay is performed on the Instrumentation Laboratory ACL TOP using the prothrombin time \(PT\) method and a factor-deficient substrate. Patient plasma is combined and incubated with a factor X-deficient substrate \(normal plasma depleted of factor X by immunoadsorption\). After a specified incubation time, a PT reagent is added to trigger the coagulation process in the mixture. Then the time to clot formation is measured optically at a wavelength of 671 nm.](#) (Owen CA Jr, Bowie EJW, Thompson JH Jr: *Diagnosis of Bleeding Disorders*. Second edition. Little, Brown and Company, Boston, MA. 1975; Meijer P, Verbruggen and Spannagi M: Chapter 33: Clotting factors and inhibitors: Assays and Interpretation. In *Laboratory Hematology Practice*. Edited by K Kottke-Marchant. Wiley Blackwell Publishing, 2012, pp 435-446)

#### PDF Report

No

#### Day(s) and Time(s) Test Performed

Monday through Friday

#### Analytic Time

1 day

#### Maximum Laboratory Time

3 days

#### Specimen Retention Time

7 days

#### Performing Laboratory Location

Rochester

---

**Fees and Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been modified from the manufacturer's instructions. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

85260

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
F_10	Coag Factor X Assay, P	3218-5

Result ID	Test Result Name	Result LOINC Value
F_10	Coag Factor X Assay, P	3218-5