

**Overview****Method Name**

High Performance Liquid Chromatography with Ultraviolet Detection (HPLC-UV)

**NY State Available**

Yes

**Specimen****Specimen Type**

Varies

**Specimen Required****Submit only 1 of the following specimens:****Plasma**

Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

**Serum**

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 2 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.3 mL

**Reject Due To**

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient	72 hours	

## Clinical and Interpretive

### Reference Values

Units: ng/mL

Plasma insulin concentrations have been shown to increase only when plasma glipizide concentrations exceeded 200 ng/mL.

Toxic range has not been established.

## Performance

### PDF Report

No

### Day(s) and Time(s) Test Performed

Monday through Sunday

### Analytic Time

7 days

### Maximum Laboratory Time

9 - 10 days

### Specimen Retention Time

2 weeks

### Performing Laboratory Location

Medtox Laboratories, Inc.

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### CPT Code Information

80299

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FGLIP	Glipizide (Glucotrol)	10539-5



---

Result ID	Test Result Name	Result LOINC Value
Z1180	Glipizide	10539-5