



# Test Definition: CBC

Complete Blood Cell Count (CBC) with  
Differential, Blood

## Overview

### Useful For

Screening tool to confirm a hematologic disorder, to establish or rule out a diagnosis, to detect an unsuspected hematologic disorder, or to monitor effects of radiation or chemotherapy

### Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
DIFFR	Morphology Eval (special Smear)	No	No
PINTP	Peripheral Smear Interpretation	No	No

### Testing Algorithm

A morphology evaluation is performed at an additional charge when the automated 5-part differential is invalid or when abnormal results are identified upon microscopic examination.

A peripheral blood smear review is performed at an additional charge when pathologist expertise is needed for clinically significant diagnosis.

### Method Name

RF/DC (Radio Frequency/Direct Current) Detection/Hydrodynamic Focusing (DC Detection), Flow Cytometry (Using a Semiconductor Laser)/Sodium Lauryl Sulfate (SLS) Hemoglobin

### NY State Available

No

## Specimen

### Specimen Type

Whole Blood EDTA

### Specimen Required

**Container/Tube:** Lavender top (EDTA)

**Specimen Volume:** 3 mL

### Forms

If not ordering electronically, complete, print, and send a [Kidney Transplant Test Request](#) with the specimen.

### Specimen Minimum Volume

1.5 mL

### Reject Due To

Gross hemolysis	Reject
Other	Clotted

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA	Refrigerated (preferred)	48 hours	
	Ambient	24 hours	

## Clinical & Interpretive

### Clinical Information

RBCs, WBCs, and platelets are produced in the bone marrow and released into the peripheral blood. The primary function of the RBC is to deliver oxygen to tissues. WBCs are key components of the immune system. Platelets play a vital role in blood clotting.

Mean corpuscular volume (MCV) is a measure of the size of the average RBC. Anemias are characterized as microcytic (MCV <80), macrocytic (MCV >100), or normocytic. The red cell distribution width (RDW) is a measure of the degree of variation in RBC size (anisocytosis). RDW may be helpful in distinguishing between some anemias. For example, iron deficiency anemia is characterized by a high RDW, while thalassemia is characterized by a low RDW.

These counts are used as clinical guides in the diagnosis or monitoring of many diseases.

### Reference Values

#### RED BLOOD CELL COUNT (RBC)

Males:

0-14 days: 4.10-5.55 x 10<sup>(12)</sup>/L

15 days-4 weeks: 3.16-4.63 x 10<sup>(12)</sup>/L

5 weeks-7 weeks: 3.02-4.22 x 10<sup>(12)</sup>/L

8 weeks-5 months: 3.43-4.80 x 10<sup>(12)</sup>/L

6 months-23 months: 4.03-5.07 x 10<sup>(12)</sup>/L

24 months-35 months: 3.89-4.97 x 10<sup>(12)</sup>/L

3-5 years: 4.00-5.10 x 10<sup>(12)</sup>/L

6-10 years: 4.10-5.20 x 10<sup>(12)</sup>/L

11-14 years: 4.20-5.30 x 10<sup>(12)</sup>/L

15-17 years: 4.30-5.70 x 10<sup>(12)</sup>/L

---

Adults: 4.35-5.65 x 10<sup>12</sup>/L

Females:

0-14 days: 4.12-5.74 x 10<sup>12</sup>/L

15 days-4 weeks: 3.32-4.80 x 10<sup>12</sup>/L

5 weeks-7 weeks: 2.93-3.87 x 10<sup>12</sup>/L

8 weeks-5 months: 3.45-4.75 x 10<sup>12</sup>/L

6 months-23 months: 3.97-5.01 x 10<sup>12</sup>/L

24 months-35 months: 3.84-4.92 x 10<sup>12</sup>/L

3-5 years: 4.00-5.10 x 10<sup>12</sup>/L

6-10 years: 4.10-5.20 x 10<sup>12</sup>/L

11-14 years: 4.10-5.10 x 10<sup>12</sup>/L

15-17 years: 3.80-5.00 x 10<sup>12</sup>/L

Adults: 3.92-5.13 x 10<sup>12</sup>/L

HEMOGLOBIN

Males:

0-14 days: 13.9-19.1 g/dL

15 days-4 weeks: 10.0-15.3 g/dL

5 weeks-7 weeks: 8.9-12.7 g/dL

8 weeks-5 months: 9.6-12.4 g/dL

6 months-23 months: 10.1-12.5 g/dL

24 months-35 months: 10.2-12.7 g/dL

3-5 years: 11.4-14.3 g/dL

6-8 years: 11.5-14.3 g/dL

9-10 years: 11.8-14.7 g/dL

11-14 years: 12.4-15.7 g/dL

15-17 years: 13.3-16.9 g/dL

Adults: 13.2-16.6 g/dL

Females:

0-14 days: 13.4-20.0 g/dL

15 days-4 weeks: 10.8-14.6 g/dL

5 weeks-7 weeks: 9.2-11.4 g/dL

8 weeks-5 months: 9.9-12.4 g/dL

6 months-35 months: 10.2-12.7 g/dL

3-5 years: 11.4-14.3 g/dL

6-8 years: 11.5-14.3 g/dL

9-10 years: 11.8-14.7 g/dL

11-17 years: 11.9-14.8 g/dL

Adults: 11.6-15.0 g/dL

HEMATOCRIT

Males:

---

0-14 days: 39.8-53.6%  
15 days-4 weeks: 30.5-45.0%  
5 weeks-7 weeks: 26.8-37.5%  
8 weeks-5 months: 28.6-37.2%  
6 months-23 months: 30.8-37.8%  
24 months-35 months: 31.0-37.7%  
3-7 years: 34-42%  
8-11 years: 35-43%  
12-15 years: 38-47%  
16-17 years: 40-50%  
Adults: 38.3-48.6%

**Females:**

0-14 days: 39.6-57.2%  
15 days-4 weeks: 32.0-44.5%  
5 weeks-7 weeks: 27.7-35.1%  
8 weeks-5 months: 29.5-37.1%  
6 months-23 months: 30.9-37.9%  
24 months-35 months: 31.2-37.8%  
3-7 years: 34-42%  
8-17 years: 35-43%  
Adults: 35.5-44.9%

**MEAN CORPUSCULAR VOLUME (MCV)****Males:**

0-14 days: 91.3-103.1 fL  
15 days-4 weeks: 89.4-99.7 fL  
5 weeks-7 weeks: 84.3-94.2 fL  
8 weeks-5 months: 74.1-87.5 fL  
6 months-23 months: 69.5-81.7 fL  
24 months-35 months: 71.3-84.0 fL  
3-5 years: 77.2-89.5 fL  
6-11 years: 77.8-91.1 fL  
12-14 years: 79.9-93.0 fL  
15-17 years: 82.5-98.0 fL  
Adults: 78.2-97.9 fL

**Females:**

0-14 days: 92.7-106.4 fL  
15 days-4 weeks: 90.1-103.0 fL  
5 weeks-7 weeks: 83.4-96.4 fL  
8 weeks-5 months: 74.8-88.3 fL  
6 months-23 months: 71.3-82.6 fL  
24 months-35 months: 72.3-85.0 fL

---

3-5 years: 77.2-89.5 fL  
6-11 years: 77.8-91.1 fL  
12-14 years: 79.9-93.0 fL  
15-17 years: 82.5-98.0 fL  
Adults: 78.2-97.93 fL

**RED CELL DISTRIBUTION WIDTH (RDW)****Males:**

0-14 days: 14.8-17.0%  
15 days-4 weeks: 14.3-16.8%  
5 weeks-7 weeks: 13.8-16.1%  
8 weeks-5 months: 12.4-15.3%  
6 months-23 months: 12.9-15.6%  
24 months-35 months: 12.5-14.9%  
3-5 years: 11.3-13.4%  
6-17 years: 11.4-13.5%  
Adults: 11.8-14.5%

**Females:**

0-14 days: 14.6-17.3%  
15 days-4 weeks: 14.4-16.2%  
5 weeks-7 weeks: 13.6-15.8%  
8 weeks-5 months: 12.2-14.3%  
6 months-23 months: 12.7-15.1%  
24 months-35 months: 12.4-14.9%  
3-5 years: 11.3-13.4%  
6-17 years: 11.4-13.5%  
Adults: 12.2-16.1%

**WHITE BLOOD CELL COUNT (WBC)****Males:**

0-14 days: 8.0-15.4 x 10<sup>(9)</sup>/L  
15 days-4 weeks: 7.8-15.9 x 10<sup>(9)</sup>/L  
5 weeks-7 weeks: 8.1-15.0 x 10<sup>(9)</sup>/L  
8 weeks-5 months: 6.5-13.3 x 10<sup>(9)</sup>/L  
6 months-23 months: 6.0-13.5 x 10<sup>(9)</sup>/L  
24 months-35 months: 5.1-13.4 x 10<sup>(9)</sup>/L  
3-5 years: 4.4-12.9 x 10<sup>(9)</sup>/L  
6-17 years: 3.8-10.4 x 10<sup>(9)</sup>/L  
Adults: 3.4-9.6 x 10<sup>(9)</sup>/L

**Females:**

0-14 days: 8.2-14.6 x 10<sup>(9)</sup>/L  
15 days-4 weeks: 8.4-14.4 x 10<sup>(9)</sup>/L

---

5 weeks-7 weeks:  $7.1-14.7 \times 10(9)/L$   
8 weeks-5 months:  $6.0-13.3 \times 10(9)/L$   
6 months-23 months:  $6.5-13.0 \times 10(9)/L$   
24 months-35 months:  $4.9-13.2 \times 10(9)/L$   
3-5 years:  $4.4-12.9 \times 10(9)/L$   
6-17 years:  $3.8-10.4 \times 10(9)/L$   
Adults:  $3.4-9.6 \times 10(9)/L$

**PLATELETS****Males:**

0-14 days:  $218-419 \times 10(9)/L$   
15 days-4 weeks:  $248-586 \times 10(9)/L$   
5 weeks-7 weeks:  $229-562 \times 10(9)/L$   
8 weeks-5 months:  $244-529 \times 10(9)/L$   
6 months-23 months:  $206-445 \times 10(9)/L$   
24 months-35 months:  $202-403 \times 10(9)/L$   
3-5 years:  $187-445 \times 10(9)/L$   
6-9 years:  $187-400 \times 10(9)/L$   
10-13 years:  $177-381 \times 10(9)/L$   
14-17 years:  $139-320 \times 10(9)/L$   
Adults:  $135-317 \times 10(9)/L$

**Females:**

0-14 days:  $144-449 \times 10(9)/L$   
15 days-4 weeks:  $279-571 \times 10(9)/L$   
5 weeks-7 weeks:  $331-597 \times 10(9)/L$   
8 weeks-5 months:  $247-580 \times 10(9)/L$   
6 months-23 months:  $214-459 \times 10(9)/L$   
24 months-35 months:  $189-394 \times 10(9)/L$   
3-5 years:  $187-445 \times 10(9)/L$   
6-9 years:  $187-400 \times 10(9)/L$   
10-13 years:  $177-381 \times 10(9)/L$   
14-17 years:  $158-362 \times 10(9)/L$   
Adults:  $157-371 \times 10(9)/L$

**NEUTROPHILS****Males:**

0-14 days:  $1.60-6.06 \times 10(9)/L$   
15 days-4 weeks:  $1.18-5.45 \times 10(9)/L$   
5 weeks-7 weeks:  $0.83-4.23 \times 10(9)/L$   
8 weeks-5 months:  $0.97-5.45 \times 10(9)/L$   
6 months-23 months:  $1.19-7.21 \times 10(9)/L$   
24 months-35 months:  $1.54-7.92 \times 10(9)/L$   
3-5 years:  $1.60-7.80 \times 10(9)/L$

---

6-16 years: 1.40-6.10 x 10<sup>9</sup>/L

17 years: 1.80-7.20 x 10<sup>9</sup>/L

Adults: 1.56-6.45 x 10<sup>9</sup>/L

Females:

0-14 days: 1.73-6.75 x 10<sup>9</sup>/L

15 days-4 weeks: 1.23-4.80 x 10<sup>9</sup>/L

5 weeks-7 weeks: 1.00-4.68 x 10<sup>9</sup>/L

8 weeks-5 months: 1.04-7.20 x 10<sup>9</sup>/L

6 months-23 months: 1.27-7.18 x 10<sup>9</sup>/L

24 months-35 months: 1.60-8.29 x 10<sup>9</sup>/L

3-5 years: 1.60-7.80 x 10<sup>9</sup>/L

6-14 years: 1.50-6.50 x 10<sup>9</sup>/L

15-17 years: 2.00-7.40 x 10<sup>9</sup>/L

Adults: 1.56-6.45 x 10<sup>9</sup>/L

LYMPHOCYTES

Males:

0-14 days: 2.07-7.53 x 10<sup>9</sup>/L

15 days-4 weeks: 2.11-8.38 x 10<sup>9</sup>/L

5 weeks-7 weeks: 2.47-7.95 x 10<sup>9</sup>/L

8 weeks-5 months: 2.45-8.89 x 10<sup>9</sup>/L

6 months-23 months: 1.56-7.83 x 10<sup>9</sup>/L

24 months-35 months: 1.13-5.52 x 10<sup>9</sup>/L

3-5 years: 1.60-5.30 x 10<sup>9</sup>/L

6-11 years: 1.40-3.90 x 10<sup>9</sup>/L

12-17 years: 1.00-3.20 x 10<sup>9</sup>/L

Adults: 0.95-3.07 x 10<sup>9</sup>/L

Females:

0-14 days: 1.75-8.00 x 10<sup>9</sup>/L

15 days-4 weeks: 2.42-8.20 x 10<sup>9</sup>/L

5 weeks-7 weeks: 2.29-9.14 x 10<sup>9</sup>/L

8 weeks-5 months: 2.14-8.99 x 10<sup>9</sup>/L

6 months-23 months: 1.52-8.09 x 10<sup>9</sup>/L

24 months-35 months: 1.25-5.77 x 10<sup>9</sup>/L

3-5 years: 1.60-5.30 x 10<sup>9</sup>/L

6-11 years: 1.40-3.90 x 10<sup>9</sup>/L

12-17 years: 1.00-3.20 x 10<sup>9</sup>/L

Adults: 0.95-3.07 x 10<sup>9</sup>/L

MONOCYTES

Males:

0-14 days: 0.52-1.77 x 10<sup>9</sup>/L

---

15 days-4 weeks:  $0.28-1.38 \times 10^9/L$   
5 weeks-7 weeks:  $0.28-1.05 \times 10^9/L$   
8 weeks-5 months:  $0.28-1.07 \times 10^9/L$   
6 months-23 months:  $0.25-1.15 \times 10^9/L$   
24 months-35 months:  $0.19-0.94 \times 10^9/L$   
3-5 years:  $0.30-0.90 \times 10^9/L$   
6-17 years:  $0.20-0.80 \times 10^9/L$   
Adults:  $0.26-0.81 \times 10^9/L$

**Females:**

0-14 days:  $0.57-1.72 \times 10^9/L$   
15 days-4 weeks:  $0.42-1.21 \times 10^9/L$   
5 weeks-7 weeks:  $0.28-1.21 \times 10^9/L$   
8 weeks-5 months:  $0.24-1.17 \times 10^9/L$   
6 months-23 months:  $0.26-1.08 \times 10^9/L$   
24 months-35 months:  $0.24-0.92 \times 10^9/L$   
3-5 years:  $0.30-0.90 \times 10^9/L$   
6-17 years:  $0.20-0.80 \times 10^9/L$   
Adults:  $0.26-0.81 \times 10^9/L$

**EOSINOPHILS****Males:**

0-14 days:  $0.12-0.66 \times 10^9/L$   
15 days-4 weeks:  $0.08-0.80 \times 10^9/L$   
5 weeks-7 weeks:  $0.05-0.57 \times 10^9/L$   
8 weeks-5 months:  $0.03-0.61 \times 10^9/L$   
6 months-23 months:  $0.02-0.82 \times 10^9/L$   
24 months-35 months:  $0.03-0.53 \times 10^9/L$   
3-11 years:  $0.00-0.50 \times 10^9/L$   
12-17 years:  $0.10-0.20 \times 10^9/L$   
Adults:  $0.03-0.48 \times 10^9/L$

**Females:**

0-14 days:  $0.09-0.64 \times 10^9/L$   
15 days-4 weeks:  $0.06-0.75 \times 10^9/L$   
5 weeks-7 weeks:  $0.04-0.63 \times 10^9/L$   
8 weeks-5 months:  $0.02-0.74 \times 10^9/L$   
6 months-23 months:  $0.02-0.58 \times 10^9/L$   
24 months-35 months:  $0.03-0.46 \times 10^9/L$   
3-11 years:  $0.00-0.50 \times 10^9/L$   
12-17 years:  $0.10-0.20 \times 10^9/L$   
Adults:  $0.03-0.48 \times 10^9/L$

**BASOPHILS**

**Males:**0-14 days: 0.02-0.11 x 10<sup>9</sup>/L15 days-7 weeks: 0.01-0.07 x 10<sup>9</sup>/L8 weeks-35 months: 0.01-0.06 x 10<sup>9</sup>/L3-17 years: 0.00-0.10 x 10<sup>9</sup>/LAdults: 0.01-0.08 x 10<sup>9</sup>/L**Females:**0-14 days: 0.02-0.07 x 10<sup>9</sup>/L15 days-4 weeks: 0.01-0.06 x 10<sup>9</sup>/L5 weeks-7 weeks: 0.01-0.05 x 10<sup>9</sup>/L8 weeks-5 months: 0.01-0.07 x 10<sup>9</sup>/L6 months-35 months: 0.01-0.06 x 10<sup>9</sup>/L3-17 years: 0.00-0.10 x 10<sup>9</sup>/LAdults: 0.01-0.08 x 10<sup>9</sup>/L**Interpretation**

Results outside of normal value ranges may reflect a primary disorder of the cell-producing organs or an underlying disease. Results should be interpreted in conjunction with the patient's clinical picture and appropriate additional testing performed.

**Cautions**

Questionable results are detected by in-house checking criteria based on quantitative and qualitative parameters.

Defined laboratory instrument or patient flagging criteria may trigger peripheral blood smear review resulting in consultant readout of hematologic abnormalities.

**Clinical Reference**

1. Adel K, Raizman J, Chen Y, et al: Complex biological profile of hematologic markers across pediatric, adult, and geriatric ages: establishment of robust pediatric and adult reference intervals on the basis of the Canadian Health Measures Survey. *Clin Chem* 2015;61:8
2. CLSI. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline, Third Edition. CLSI document EP28-A3c. Wayne, PA, Clinical and Laboratory Standards Institute, 2008
3. Soldin J, Brugnara C, Wong EC: Pediatric Reference Intervals. Fifth Edition. AACC Press. Washington, DC, 2005. ISBN 1-594250-32-4

**Performance****Method Description**

The Sysmex XE-5000 performs hematology analysis according to the radio frequency/direct current (RF/DC) detection method. Hydrodynamic focusing (DC detection), flow cytometry method (using a semiconductor laser), and sodium lauryl sulfate (SLS)-hemoglobin method. (Instruction manual: Automated Hematology Analyzer XE-5000 [North American Edition], Code No. 461-2638-6, rev: June 2008, Copyright 2006-2008 by SYSMEX CORPORATION)

**PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

Same day/1 day

**Specimen Retention Time**

1 week

**Performing Laboratory Location**

Mayo Clinic Health System in Mankato

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been cleared, approved, or is exempt by the US Food and Drug Administration and is used per manufacturer's instructions. Performance characteristics were verified by Mayo Clinic in a manner consistent with CLIA requirements.

**CPT Code Information**

85025

85007 (if appropriate)

85060 (if appropriate)

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
CBC	CBC with Differential, B	57021-8

Result ID	Test Result Name	Result LOINC® Value
HGB	Hemoglobin	718-7
HCT	Hematocrit	4544-3
RBC	Erythrocytes	789-8
MCV	MCV	787-2

## Test Definition: CBC

Complete Blood Cell Count (CBC) with  
Differential, Blood

RDW	RDW	788-0
WBC	Leukocytes	6690-2
NEUAA	Neutrophils, Abs	751-8
LYMAA	Lymphocytes, Abs	731-0
MONAA	Monocytes, Abs	742-7
EOSAA	Eosinophils, Abs	711-2
BASAA	Basophils, Abs	704-7
PLTC	Platelet Count	777-3