

Overview

Useful For

Detecting circulating immune complexes

Method Name

Semi-Quantitative Enzyme-Linked Immunosorbent Assay

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Collection Container/Tube: Red top or serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL serum

Collection Instructions:

1. Allow complete clotting of red blood cells (up to 1 hour).
2. Within 30 minutes after clotting is complete, centrifuge and aliquot serum into a plastic vial.
3. Freeze immediately.
4. Ship frozen.

Note: Critical frozen. Separate specimens must be submitted when multiple tests are ordered.

Specimen Minimum Volume

Serum: 0.3 mL

Reject Due To

Gross hemolysis:	Reject
Thawing:	Cold reject; Warm reject
Non-frozen specimen	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Serum	Frozen	14 days	
-------	--------	---------	--

Clinical & Interpretive

Reference Values

Less than or equal to 3.9 ug Eq/mL

Interpretation

Less than or equal to 3.9 ug Eq/mL is considered negative for circulating complement binding immune complexes.

Circulating immune complexes may be found without any evident pathology and positive results do not necessarily implicate the immune complex in a disease process.

Performance

PDF Report

No

Day(s) Performed

Saturday, Monday, and Thursday

Report Available

4 to 15 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86332

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCQBA	CIC, C1q Binding	27831-7

Result ID	Test Result Name	Result LOINC® Value
Z2636	CIC, C1q Binding	27831-7