



# Test Definition: FFPG

Fruit Panel IgG

## Overview

### Method Name

Enzyme immunoassay (FEIA)

### NY State Available

No

## Specimen

### Specimen Type

Serum

### Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1.5 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

1.5 mL

### Reject Due To

|           |                       |
|-----------|-----------------------|
| Hemolysis | Mild OK; Gross reject |
| Lipemia   | Mild OK; Gross reject |
| Icterus   | Mild OK; Gross reject |
| Other     | NA                    |

### Specimen Stability Information

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Ambient                  | 7 days   |                   |
|               | Frozen                   | 365 days |                   |

## Clinical & Interpretive

### Reference Values

Reference ranges have not been established for food-specific IgG tests. The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet

elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

3 to 7 days

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

### CPT Code Information

86001 x16

### LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------|--------------------|
| FFPG    | Fruit Panel IgG | Not Provided       |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| Z3968     | Apple IgG        | 60349-8             |
| Z3969     | Apricot IgG      | 60350-6             |
| Z3970     | Banana IgG       | 60352-2             |
| Z3971     | Blueberry IgG    | 58035-7             |

|       |                |         |
|-------|----------------|---------|
| Z3972 | Cranberry IgG  | 63104-4 |
| Z3973 | Grape IgG      | 60376-1 |
| Z3974 | Grapefruit IgG | 60377-9 |
| Z3975 | Lemon IgG      | 60381-1 |
| Z3976 | Orange IgG     | 51904-1 |
| Z3977 | Papaya IgG     | 63131-7 |
| Z3978 | Peach IgG      | 60390-2 |
| Z3979 | Pear IgG       | 60391-0 |
| Z3980 | Pineapple IgG  | 60394-4 |
| Z3981 | Plum IgG       | 60395-1 |
| Z3982 | Raspberry IgG  | 63083-0 |
| Z3983 | Strawberry IgG | 40835-1 |