

# **Test Definition: FSFM**

Soluble Fibrin Monomer

# **Overview**

#### **Method Name**

Qualitative Hemagglutination

#### **NY State Available**

Yes

# **Specimen**

# **Specimen Type**

Plasma Na Cit

# **Specimen Required**

Collection Container/Tube: Light-blue top (3.2% sodium citrate)

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL Platelet-poor plasma

**Collection Instructions:** 

- 1. Centrifuge at 1700 x g for 15 minutes.
- 2. Aliquot 1.5 platelet-poor plasma into plastic vial making sure not to disturb the buffy coat.
- 3. Send frozen.

Note: Critical frozen. Separate specimens must be submitted when multiple tests are ordered.

# **Specimen Minimum Volume**

Platelet-poor plasma: 1 mL

# **Reject Due To**

Hemolysis	Reject
Clotted	Reject
Glass vials	Reject
Serum, EDTA	Reject
plasma	

# **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

# **Clinical & Interpretive**



# **Test Definition: FSFM**

Soluble Fibrin Monomer

#### **Clinical Information**

Refer to https://ltd.aruplab.com/

#### **Reference Values**

Negative

#### **Cautions**

Activation of specimen due to poor venipuncture technique may cause false-positive results.

#### **Performance**

# **PDF Report**

No

# Day(s) Performed

Monday through Sunday

#### Report Available

1 to 6 days

# **Performing Laboratory Location**

**ARUP Laboratories** 

# **Fees & Codes**

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

#### **CPT Code Information**

85366

# **LOINC®** Information

FSFM Soluble Fibrin Monomer 40702-3	Test ID	Test Order Name	Order LOINC® Value
	FSFM	Soluble Fibrin Monomer	40702-3

Result ID	Test Result Name	Result LOINC® Value
FSFM	Soluble Fibrin Monomer	40702-3