

**Test Definition: SP7CC** 

Septin-7 Antibody, Cell-Binding Assay, Spinal

Fluid

#### Overview

#### **Useful For**

Detecting septin-7 IgG by cell-binding assay using cerebrospinal fluid specimens

# **Testing Algorithm**

If the indirect immunofluorescence pattern suggests septin-7, then this test and septin-7 antibody IFA titer will be performed at an additional charge.

#### **Method Name**

Only orderable as a reflex. For more information see:

- -ENC2 / Encephalopathy, Autoimmune /Paraneoplastic Evaluation, Spinal Fluid
- -MDC2 / Movement Disorder, Autoimmune /Paraneoplastic Evaluation, Spinal Fluid
- -MAC1 / Myelopathy, Autoimmune/Paraneoplastic Evaluation, Spinal Fluid

Cell-Binding Assay (CBA)

#### **NY State Available**

Yes

## Specimen

## **Specimen Type**

**CSF** 

# **Specimen Required**

Only orderable as a reflex. For more information see:

- -ENC2 / Encephalopathy, Autoimmune /Paraneoplastic Evaluation, Spinal Fluid
- -MDC2 / Movement Disorder, Autoimmune /Paraneoplastic Evaluation, Spinal Fluid
- -MAC1 / Myelopathy, Autoimmune/Paraneoplastic Evaluation, Spinal Fluid

Container/Tube: Sterile vial Specimen Volume: 1.5 mL

#### Specimen Minimum Volume

See Specimen Required

# Reject Due To

Gross	Reject
hemolysis	



# **Test Definition: SP7CC**

Septin-7 Antibody, Cell-Binding Assay, Spinal

Fluid

Gross lipemia	Reject
Gross icterus	Reject

# **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
CSF	Refrigerated (preferred)	28 days	
	Ambient	72 hours	
	Frozen	28 days	

# **Clinical & Interpretive**

#### **Clinical Information**

Neurological phenotypes for septin-7 IgG positive patients include encephalopathy, myelopathy, encephalomyelopathy, painful myelopolyradiculopathy, and episodic ataxia. Psychiatric symptoms are also common with encephalopathic symptoms. Septin-7 IgG is also associated with cancer. Positive response to immunotherapy.

#### **Reference Values**

Only orderable as a reflex. For more information see:

- -ENC2 / Encephalopathy, Autoimmune /Paraneoplastic Evaluation, Spinal Fluid
- -MDC2 / Movement Disorder, Autoimmune /Paraneoplastic Evaluation, Spinal Fluid
- -MAC1 / Myelopathy, Autoimmune/Paraneoplastic Evaluation, Spinal Fluid

Negative

#### Interpretation

Seropositivity for septin antibodies by indirect immunofluorescence assay is consistent with a diagnosis of autoimmune disease of the central nervous system. Cell-binding assay (CBA) testing for septin-7 IgG is required to confirm the diagnosis. Seropositivity for Septin-7 IgG by CBA confirms a diagnosis of autoimmune disease of the central nervous system. A paraneoplastic cause should be considered in patients with septin-7-IgG.

### **Cautions**

Negative results for septin-7 IgG by cell-binding assay do not exclude neurological autoimmunity or cancer.

## **Clinical Reference**

Honorat JA, Miske R, Scharf M, et al: 416. Neuronal septin autoimmunity: Differentiated serological profiles and clinical findings. Ann Neurol. 2020 Oct;88(Suppl 25):S55. Abstract

# **Performance**



# **Test Definition: SP7CC**

Septin-7 Antibody, Cell-Binding Assay, Spinal

Fluid

## **Method Description**

Patient specimen is applied to a composite slide containing transfected and nontransfected HEK-293 cells. After incubation and washing, fluorescein-conjugated goat-antihuman IgG is applied to detect the presence of patient IgG binding.(Package insert: IIFT: Neurology Mosaics, Instructions for the indirect immunofluorescence test. EUROIMMUN; FA\_112d-1\_A\_UK\_C13, 02/2019)

## **PDF Report**

No

# Day(s) Performed

Monday through Sunday

# **Report Available**

5 to 10 days

# **Specimen Retention Time**

28 days

# **Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Main Campus

#### Fees & Codes

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

# **Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

## **CPT Code Information**

86255

# **LOINC®** Information

	NC® Value
SP7CC Septin-7 CBA, CSF 101466-1	

Result ID	Test Result Name	Result LOINC® Value
615872	Septin-7 CBA, CSF	101466-1