



# Test Definition: CALFM

CALR Mutation Analysis, Myeloproliferative Neoplasm (MPN), Bone Marrow

## Overview

### Useful For

[Diagnosis or differential diagnosis of myeloproliferative disorders by CALR gene sequencing using bone marrow specimens](#)

### Method Name

Polymerase Chain Reaction (PCR) and Fragment Analysis

### NY State Available

No

## Specimen

### Specimen Type

Bone Marrow

### Shipping Instructions

Specimen must arrive within 7 days of collection.

### Necessary Information

The following information is required:

1. Pertinent clinical history
2. Clinical or morphologic suspicion
3. Date of collection
4. Specimen source

### Specimen Required

**Container/Tube:** Lavender top (EDTA)

**Specimen Volume:** 3 mL

#### Collection Instructions:

1. Invert several times to mix bone marrow.
2. Send bone marrow specimen in original tube. **Do not aliquot.**
3. Label specimen as bone marrow.

### Specimen Minimum Volume

1 mL

### Reject Due To

Gross	Reject
-------	--------

hemolysis	
Moderately to severely clotted	Reject

## Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Bone Marrow	Ambient (preferred)	7 days	
	Refrigerated	7 days	

## Clinical & Interpretive

### Clinical Information

Mutations in the *JAK2*, *CALR* and *MPL* genes are considered driver events in the BCR-ABL1 negative myeloproliferative neoplasms (MPN) including polycythemia vera (PV), primary myelofibrosis (PMF) and essential thrombocythemia (ET). The *JAK2* V617F mutation occurs in 95% to 98% of patients with PV, 50% to 60% of patients with PMF and 50% to 60% of patients with ET respectively at diagnosis. Other *JAK2* mutations in exon 12 to 15 occur in the remaining patients with PV. Mutations in the *CALR* gene occur in 20% to 30% of patients with PMF and 20% to 30% of patients with ET at diagnosis. A 52 base pairs (bp) deletion (53%) and a 5 bp deletion (32%) are the most common mutations in the *CALR* gene while other types of mutations may occur in the remaining cases. *MPL* exon 10 mutations occur in 5% to 10% of patients with PMF and 5% to 10% of patients with ET. Mutations in *JAK2*, *CALR* and *MPL* are mutually exclusive. The *JAK2* V617F mutation is detected by quantitative polymerase chain reaction (qPCR). The *CALR* mutations are detected by PCR targeting the exon 9. The *MPL* mutations in exon 10 are detected by Sanger sequencing. All mutations in *JAK2*, *CALR* and *MPL* can also be detected by next generation sequencing (NGS). In addition to the mutations in *JAK2*, *CALR* and *MPL*, mutations in many other genes including *ASXL1*, *TET2*, *DNMT3A*, *SRSF2*, *SF3B1*, *U2AF1*, *ZRSR2*, *EZH2*, *IDH1*, *IDH2*, *CBL*, *KRAS*, *NRAS*, *STAG2*, and *TP53* can occur in MPN. These additional mutations are more frequent in PMF and advanced disease, as compared to PV and ET. It is known that mutations in the *ASXL1*, *SRSF2*, *U2AF1*, *EZH2*, *IDH1* and *IDH2* are correlated with a poor prognostic risk. While a single gene test on *JAK2*, *CALR* and *MPL* can be clinically useful, all above mentioned gene mutations can be detected by NGS.

### Reference Values

An interpretive report will be provided.

### Interpretation

The results will be reported as 1 of the 3 states if DNA amplification is successful (see Cautions):

- Positive. A deletion-insertion-type mutation was detected in *CALR*, exon 9.
- Negative. No deletion or insertion was detected in *CALR*, exon 9.
- Equivocal. A small amplicon suspicious for a deletion-insertion type mutation was detected in *CALR*, exon 9.

Positive mutation status is highly suggestive of a myeloid neoplasm but must be correlated with clinical and other laboratory and morphologic features for definitive diagnosis.

---

Negative mutation status does not exclude the presence of a myeloproliferative neoplasm or other neoplastic disorders.

**Cautions**

A positive result is not specific for a particular myeloproliferative neoplasm (MPN) diagnosis, and clinicopathologic correlation is necessary in all cases.

A negative result does not exclude the presence of MPN or other neoplastic processes.

This test is a fragment analysis assay and only detects deletions-insertions (delins). It will not detect point mutations. However, all reported disease-causing mutations in MPN described to date are insertions and/or deletions.

This test may not differentiate between out-of-frame and in-frame delins in rare cases. However, in-frame delin mutations are very rare (<0.5%) and have only been reported in few healthy individuals and myeloproliferative neoplasm patients with *JAK2V617F* mutation or out-of-frame *CALR* mutation. Most of the rare in-frame delins are considered germline variants and represent benign alterations (ie, polymorphisms).

Infrequently, amplification failure can be encountered in a sample, due to inadequate DNA, poor DNA quality, or a polymerase chain reaction inhibitor. In these circumstances, the assay will be reattempted and if persistently unsuccessful, the report will be issued with an "Invalid" result.

**Clinical Reference**

1. Klampfl T, Gisslinger H, Harutyunyan AS, et al. Somatic mutations of calreticulin in myeloproliferative neoplasms. *N Engl J Med*. 2013;369(25):2379-2390
2. Rumi E, Pietra D, Ferretti V, et al. JAK2 or CALR mutation status defines subtypes of essential thrombocythemia with substantially different clinical course and outcomes. *Blood*. 2014;123(10):1544-1551
3. Greenfield G, McMullin MF, Mills K. Molecular pathogenesis of the myeloproliferative neoplasms. *J Hematol Oncol*. 2021;14(1):103
4. Khoury JD, Solary E, Abla O, et al. The 5th edition of the World Health Organization Classification of Haematolymphoid Tumours: Myeloid and Histiocytic/Dendritic Neoplasms. *Leukemia*. 2022;36(7):1703-1719.  
doi:10.1038/s41375-022-01613-1

**Performance****Method Description**

Polymerase chain reaction (PCR) amplification of *CALR* exon 9 is performed on DNA isolated from the patient sample. The PCR product is then run on an ABI 3130xl Genetic Analyzer for fragment analysis to detect insertions and deletions. An unmutated *CALR* will show an amplicon at 266 base pairs (bp), a mutated *CALR* with insertion will show an amplicon greater than 266 bp, and a mutated *CALR* with deletion will show an amplicon smaller than 266 bp. This assay has an analytical sensitivity of approximately 6% (ie, 6 mutation-containing cells in 100 total cells) in most mutation types, except for the rare type of 1-bp deletion, which has a sensitivity of approximately 20%. (Unpublished Mayo method)

**PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

1 to 8 days

**Specimen Retention Time**

Bone marrow: 2 weeks; Extracted DNA: 1 year

**Performing Laboratory Location**

Mayo Clinic Jacksonville Clinical Lab

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

**CPT Code Information**

81219

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
CALFM	CALR Gene Mutation, Exon 9, BM	77174-1

Result ID	Test Result Name	Result LOINC® Value
618655	Final Diagnosis	69047-9
618656	Method	85069-3
618657	Signing Pathologist	18771-6
618659	Additional Information	48767-8
618660	Disclaimer	62364-5