

MPL Exon 10 Mutation Detection, Bone
Marrow

Overview

Useful For

Diagnosis or differential diagnosis of myeloproliferative disorders by MPL gene analysis using bone marrow specimens

Method Name

Sanger Sequencing

NY State Available

No

Specimen

Specimen Type

Bone Marrow

Shipping Instructions

Specimen must arrive within 7 days of collection.

Necessary Information

The following information is required:

- 1. Pertinent clinical history
- 2. Clinical or morphologic suspicion
- 3. Date of collection
- 4. Specimen source

Specimen Required

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL Collection Instructions:

- 1. Invert several times to mix bone marrow.
- 2. Send bone marrow specimen in original tube. Do not aliquot.
- 3. Label specimen as bone marrow.

Specimen Minimum Volume

1 mL

Reject Due To

Gross	Reject
hemolysis	



MPL Exon 10 Mutation Detection, Bone
Marrow

Moderately to	Reject
severely	
clotted	

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Bone Marrow	Ambient (preferred)	7 days	
	Refrigerated	7 days	

Clinical & Interpretive

Clinical Information

Mutations in the JAK2, CALR and MPL genes are considered driver events in the BCR::ABL1 negative myeloproliferative neoplasms (MPN) including polycythemia vera (PV), primary myelofibrosis (PMF) and essential thrombocythemia (ET). The JAK2 V617F mutation occurs in 95% to 98% of patients with PV, 50% to 60% of patients with PMF and 50% to 60% of patients with ET respectively at diagnosis. Other JAK2 mutations in exons 12 to 15 occur in the remaining patients with PV. Mutations in the CALR gene occur in 20% to 30% of patients with PMF and 20% to 30% of patients with ET at diagnosis. A 52 base pairs deletion (53%) and a 5 bp deletion (32%) are the most common mutations in the CALR gene while other types of mutations may occur in the remaining cases. MPL exon 10 mutations occur in 5% to 10% of patients with PMF and 5% to 10% of patients with ET. Mutations in JAK2, CALR and MPL are mutually exclusive. The JAK2 V617F mutation is detected by quantitative polymerase chain reaction (PCR). The CALR mutations are detected by PCR targeting exon 9. The MPL mutations in exon 10 are detected by Sanger sequencing. All mutations in JAK2, CALR and MPL can also be detected by next generation of sequencing (NGS). In addition to the mutations in JAK2, CALR and MPL, mutations in many other genes including ASXL1, TET2, DNMT3A, SRSF2, SF3B1, U2AF1, ZRSR2, EZH2, IDH1, IDH2, CBL, KRAS, NRAS, STAG2, and TP53 can occur in MPN. These additional mutations are more frequent in PMF and advanced disease, as compared to PV and ET. It is known that mutations in the ASXL1, SRSF2, U2AF1, EZH2, IDH1 and IDH2 are correlated with a poor prognostic risk. While a single gene test on JAK2, CALR and MPL can be clinically useful, all above mentioned gene mutations can be detected by NGS.

Reference Values

An interpretive report will be provided.

Interpretation

The results will be reported as 1 of 2 states:

- -Negative for myeloproliferative leukemia virus oncogene (MPL) exon 10 mutation
- -Positive for MPL exon 10 mutation

If the result is positive, a description of the mutation at the nucleotide level and the altered protein sequence is reported.

Positive mutation status is highly suggestive of a myeloproliferative neoplasm but must be correlated with clinical and other laboratory features for a definitive diagnosis. Negative mutation status does not exclude the presence of a



MPL Exon 10 Mutation Detection, Bone
Marrow

myeloproliferative or other neoplasm.

Cautions

A positive result is not specific for a particular diagnosis, and clinicopathologic correlation is necessary in all cases.

A negative result does not exclude the presence of a myeloproliferative or other neoplasm.

Supportive Data

Analytical sensitivity is approximately 20% meaning there must be about 20% of the altered DNA in the specimen for reliable detection.

Clinical Reference

- 1. Klampfl T, Gisslinger H, Harutyunyan AS, et al. Somatic mutations of calreticulin in myeloproliferative neoplasms. N Engl J Med. 2013;369(25):2379-2390
- 2. Rumi E, Pietra D, Ferretti V, et al. JAK2 or CALR mutation status defines subtypes of essential thrombocythemia with substantially different clinical course and outcomes. Blood. 2014;123(10):1544-1551
- 3. Greenfield G, McMullin MF, Mills K. Molecular pathogenesis of the myeloproliferative neoplasms. J Hematol Oncol. 2021;14(1):103
- 4. Khoury JD, Solary E, Abla O, et al. The 5th edition of the World Health Organization classification of haematolymphoid tumors: myeloid and histiocytic/dendritic neoplasms. Leukemia 2022; 36:1703-1719.

Performance

Method Description

Genomic DNA is extracted from bone marrow, and the *MPL* exon 10 amplified using standard polymerase chain reaction. The entire exon 10 sequence is obtained using Sanger sequencing (BigDye terminator V3.1 cycle sequencing kid from Applied Bioscience) with analysis on an automated genetic analyzer. (Unpublished Mayo method)

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

1 to 8 days

Specimen Retention Time

Bone marrow: 2 weeks; Extracted DNA: 1 year

Performing Laboratory Location

Mayo Clinic Jacksonville Clinical Lab



MPL Exon 10 Mutation Detection, Bone Marrow

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

81339

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
MPLFM	MPL Exon 10 Mutation, BM	75033-1

Result ID	Test Result Name	Result LOINC® Value
618662	Final Diagnosis	69047-9
618663	Method	85069-3
618664	Signing Pathologist	18771-6
618666	Additional Information	48767-8
618667	Disclaimer	62364-5