



# Test Definition: ETGC

Ethyl Glucuronide Confirmation, Random,  
Urine

## Overview

### Useful For

Monitoring abstinence in clinical and justice system settings

### Method Name

Liquid Chromatography Tandem Mass Spectrometry (LC-MS/MS)

### NY State Available

Yes

## Specimen

### Specimen Type

Urine

### Ordering Guidance

For situations where chain of custody is required, a Chain of Custody Kit (T282) is available. For chain-of-custody testing, order ETGX / Ethyl Glucuronide Confirmation, Chain of Custody, Random, Urine.

Additional drug panels and specific requests are available: call 800-533-1710.

### Additional Testing Requirements

If urine creatinine is required or adulteration of the specimen is suspected, the following test should also be ordered, ADULT / Adulterants Survey, Random, Urine.

### Specimen Required

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Collection Container/Tube:** Plastic urine container

**Submission Container/Tube:** Plastic, 5-mL tube

**Specimen Volume:** 1 mL

#### Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.

#### Additional Information:

1. No specimen substitutions.
2. STAT requests are **not accepted** for this test.
3. Submitting less than 1 mL will compromise our ability to perform all necessary testing.

## Forms

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

[-Kidney Transplant Test Request](#)

[-Therapeutics Test Request](#) (T831)

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

**Clinical & Interpretive****Clinical Information**

Ethyl glucuronide (EtG) and ethyl sulfate (EtS) are minor metabolites of ethanol that are detectable in body fluids following alcohol consumption and, less commonly, following extraneous exposure. EtG and EtS can be detected up to 5 days in urine using a cutoff of 500 ng/mL.(1)

**Reference Values**

Negative

Cutoff concentrations by liquid chromatography tandem mass spectrometry:

Ethyl glucuronide: 250 ng/mL

Ethyl sulfate: 100 ng/mL

**Interpretation**

A positive interpretation will be given if either the ethyl glucuronide (EtG) result is greater than or equal to 250 ng/mL or the ethyl sulfate (EtS) is greater than or equal to 100 ng/mL.

A "high" positive (ie, >1000 ng/mL) may indicate:

- Heavy drinking on the same day or previously (ie, previous day or 2).
- Light drinking the same day

A "low" positive (ie, 500-1000 ng/mL) may indicate:

- Previous heavy drinking (ie, previous 1-3 days)
- Recent light drinking (ie, past 24 hours)
- Recent intense "extraneous" exposure (ie, within 24 hours or less)

A "very low" positive (ie, 100-500 ng/mL) may indicate:

- Previous heavy drinking (ie, 1-3 days)
- Previous light drinking (ie, 12-36 hours)
- Recent "extraneous" exposure(2)

**Cautions**

Incidental exposure to alcohol in many daily use products (ie, hand sanitizers, mouthwash) may result in detectable levels of ethyl glucuronide (EtG) or ethyl sulfate (EtS).

In addition, upper respiratory infections as well as beta-glucuronidase hydrolysis may lower levels of EtG but do not seem to affect EtS.(2)

EtG/EtS results should be interpreted in the context of all available clinical and behavioral information.

**Clinical Reference**

1. Reisfield GM, Goldberger BA, Crews BO, et al. Ethyl glucuronide, ethyl sulfate, and ethanol in urine after sustained exposure to an ethanol-based hand sanitizer. *J Anal Toxicol.* 2011;35(2):85-91. doi:10.1093/anatox/35.2.85
2. Substance Abuse and Mental Health Services Administration (SAMSHA) Advisory: The role of biomarkers in the treatment of alcohol use disorders, 2012 Revision. HHS; 2012;11(2):1-7. doi:10.1037/e558582006-001
3. Langman LJ, Bechtel LK, Holstege CP. Clinical toxicology. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, eds. *Tietz Textbook of Laboratory Medicine.* 7th ed. Elsevier; 2023:chap 43

**Performance****Method Description**

The urine sample is diluted with internal standard in 0.1% formic acid for detection by tandem mass spectrometry.(Unpublished Mayo method)

**PDF Report**

No

**Day(s) Performed**

Monday, Wednesday, Friday

**Report Available**

3 to 5 days

**Specimen Retention Time**

14 days

**Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Superior Drive

---

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

**CPT Code Information**

G0480

80321 (if appropriate for select payers)

[Clinical Toxicology CPT Code Client Guidance](#)**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
ETGC	Ethyl Glucuronide Confirmation, U	93705-2

Result ID	Test Result Name	Result LOINC® Value
63421	Ethyl Glucuronide	58378-1
36848	Ethyl Sulfate	58425-0
36849	Ethyl Gluc/Sulfate Interpretation	69050-3