



Test Definition: MUSK

Muscle-Specific Kinase (MuSK) Autoantibody,
Serum

Overview

Useful For

Diagnosis of autoimmune muscle-specific kinase (MuSK) myasthenia gravis

Second-order test to aid in the diagnosis of autoimmune myasthenia gravis when first-line serologic tests are negative

Establishing a quantitative baseline value for MuSK antibodies that allows comparison with future levels if weakness is worsening

Method Name

Radioimmunoassay (RIA)

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube:

Preferred: Red top

Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

Forms

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

[-Neurology Specialty Testing Client Test Request \(T732\)](#)

[-General Request \(T239\)](#)

Specimen Minimum Volume

1 mL

Reject Due To

Gross hemolysis	Reject
-----------------	--------

Gross lipemia	Reject
Gross icterus	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	72 hours	
	Frozen	28 days	

Clinical & Interpretive

Clinical Information

Fatigable weakness due to impaired synaptic transmission at the neuromuscular junction is characteristic of myasthenia gravis (MG). The diagnosis is made by clinical and electromyographic criteria. Positive autoimmune serology must be interpreted in the clinical and electrophysiological context and response to anticholinesterase medication. Most cases are autoimmune and are caused by IgG autoantibodies binding to critical postsynaptic membrane molecules (nicotinic acetylcholine receptor or its interacting proteins).(1) Autoantibody detection frequency is lowest in patients with weakness confined to extraocular muscles (71% muscle acetylcholine receptor: AChR binding).(2) Mayo Clinic Laboratories first-line serological evaluation detects muscle AChR antibody in 92% of nonimmunosuppressed patients with generalized weakness due to MG. Muscle-specific kinase (MuSK) antibody is detectable in more than one-third of those seronegative for muscle AChR antibody (less than 4% of all patients).(3) Physiologically, MuSK is involved in integrating and stabilizing AChR clusters in the motor endplate. MuSK is activated when the nerve-derived proteoglycan agrin binds to its receptor, lipoprotein-related protein 4 (LRP4). Antibodies to LRP4 itself have been described in rare patients.(1)

Six percent of nonimmunosuppressed patients with generalized MG lack demonstrable AChR or MuSK antibodies (double seronegative). Other rare autoantibodies no doubt remain to be discovered in such cases. However, as in autoimmune AChR MG and MuSK MG, testing for common organ-specific and nonorgan-specific autoantibodies is a valuable ancillary investigation in evaluating seronegative acquired generalized MG. General serological testing, coupled with family or personal history, will disclose autoimmune phenomena in 77% of those cases.(3) These disorders may include thyroid disease, type 1 diabetes, vitiligo, premature greying, rheumatoid arthritis, or lupus. Testing may also reveal antinuclear antibodies, glutamic acid decarboxylase (GAD65) antibodies, thyroperoxidase/thyroglobulin antibodies, or gastric parietal cell antibodies.(3) Objective improvement in strength following a therapeutic trial of plasmapheresis or intravenous immune globulin would justify consideration of long-term immunosuppression.

Female patients are generally affected by autoimmune MuSK MG more often than male patients. Onset can occur at any age (pediatric to older adults). Patients may derive limited benefit from anticholinesterase medication. The thymus is normal, and patients are generally not benefited by thymectomy. Antibody-lowering therapies are effective. Bulbar, facial, and respiratory weakness are prominent, and crises are common.(1,4)

Reference Values

< or =0.02 nmol/L

Interpretation

A positive result, in the appropriate clinical context, confirms the diagnosis of autoimmune muscle-specific kinase myasthenia gravis.

Seropositivity justifies consideration of immunotherapy.

Cautions

Immunosuppressant therapy is a common cause of false-seronegativity. It is, therefore, important to perform a comprehensive serological evaluation before initiating immunosuppressant therapy.

Interpretation of a patient's serological and clinical status is further complicated when characteristic signs of myasthenia gravis are obscured by a superimposed steroid-induced myopathy.

Clinical Reference

1. Li Y, Arora Y, Levin K. Myasthenia gravis: newer therapies offer sustained improvement. *Cleve Clin J Med.* 2013;80(11):711-721
2. Lennon VA: Serological profile of myasthenia gravis and distinction from the Lambert-Eaton myasthenic syndrome. *Neurology* 1997;48 (Suppl 5):S23-S27
3. Chan KH, Lachance DH, Harper CM, Lennon VA. Frequency of seronegativity in adult-acquired generalized myasthenia gravis. *Muscle Nerve.* 2007;36(5):651-658
4. Skjei KL, Lennon VA, Kuntz NL. Muscle specific kinase autoimmune myasthenia gravis in children: A case series. *Neuromuscul Disord.* 2013;23(11):874-882

Performance**Method Description**

(125)I-labeled recombinant human muscle-specific kinase (MuSK) is incubated with patient sample. Anti-human IgG is then added to form an immunoprecipitate. After washing the immunoprecipitate, the amount of (125) I-labeled antigen in the immunoprecipitate is measured using a gamma-counter. The amount of gamma emission in the precipitate is proportional to the amount of MuSK-IgG in the sample. Results are reported as units of precipitated antigen (nMol) per L of patient sample. (Lavrnic D, Losen M, Vujic A, et al. The features of myasthenia gravis with autoantibodies to MuSK. *J Neurol Neurosurg Psychiatry.* 2005;76[8]:1099-1102)

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 10 days

Specimen Retention Time

2 days

Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Main Campus

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed using an analyte specific reagent. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

86366

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
MUSK	MuSK Autoantibody, S	51716-9

Result ID	Test Result Name	Result LOINC® Value
64277	MuSK Autoantibody, S	51716-9