



Test Definition: FCBDS

Cannabidiol, Serum

Overview

Method Name

High Performance Liquid Chromatography/ Tandem Mass Spectrometry (LC-MS/MS)

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Collection Container/Tube: Red top (Serum gel/SST are not acceptable)

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL Serum

Collections Instructions:

1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Stability Information: Ambient 14 days/Refrigerated 30 days/Frozen (-20 degrees C) 2 years

Specimen Minimum Volume

Serum: 0.5 mL

Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|---------|-------------------|
| Serum Red | Refrigerated (preferred) | 30 days | |
| | Ambient | 14 days | |
| | Frozen | | |

Clinical & Interpretive

Clinical Information

Refer to www.nmslabs.com/test-catalog.

Reference Values

Reporting limit determined each analysis

Interpretation

Mean peak CBD plasma concentrations at 3-4 hours post-dose with Sativex at a low dose (5.4 mg of Delta-9 THC and 5.0 mg of Cannabidiol) were 1.6 +/- 0.4 ng/mL and at a high dose (16 mg of Delta-9 THC and 15 mg of Cannabidiol) were 6.7 +/- 2.0 ng/mL.

Following high dose 400 and 800 mg oral synthetic CBD in corn oil administration, mean peak CBD plasma concentrations occurred within 1.5-3 hours post-dose and were 181.2 +/-39.8 and 221.1 +/-35.6 ng/mL, respectively. The ratio of blood to serum or plasma concentration is unknown for this analyte.

Cautions

Substance(s) known to interfere with the identity and/or quantity of the reported result: oleamide

Performance**PDF Report**

No

Day(s) Performed

Monday through Sunday

Report Available

11 to 13 days

Performing Laboratory Location

NMS Labs

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80299

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|--------------------|--------------------|
| FCBDS | Cannabidiol, Serum | 82964-8 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| Z5547 | Cannabidiol | 82964-8 |