



# Test Definition: FFWNC

Walnut Component rJug r 3

## Overview

### Method Name

ImmunoCAP FEIA

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

**Specimen Type:** Serum

**Container/Tube:** Red or SST

**Specimen Volume:** 0.5 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

0.5 mL

### Reject Due To

|                  |                  |
|------------------|------------------|
| Gross hemolysis: | Reject           |
| Thawing:         | Warm OK; Cold OK |
| Gross lipemia:   | Reject           |
| Gross icterus:   | NA               |
| Other:           | NA               |

### Specimen Stability Information

| Specimen Type | Temperature              | Time    | Special Container |
|---------------|--------------------------|---------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days |                   |
|               | Ambient                  | 28 days |                   |
|               | Frozen                   | 84 days |                   |

## Clinical & Interpretive

**Reference Values**

&lt;0.10 kU/L

**Interpretation**

| Class | IgE (kU/L)    | Comment              |
|-------|---------------|----------------------|
| 0     | <0.10         | Negative             |
| 0/1   | 0.10 - 0.34   | Equivocal/Borderline |
| 1     | 0.35 - 0.69   | Low Positive         |
| 2     | 0.70 - 3.49   | Moderate Positive    |
| 3     | 3.50 - 17.49  | High Positive        |
| 4     | 17.50 - 49.99 | Very High Positive   |
| 5     | 50.00 - 99.99 | Very High Positive   |
| 6     | >99.99        | Very High Positive   |

**Performance****Method Description**

The ImmunoCAP FEIA method uses as the solid phase a flexible, hydrophobic cellulosic polymer to which allergen has been covalently linked. The advantage of this system is that it has a very high antigen binding capacity when compared to other systems and it has minimal non-specific binding with high total IgE.

**PDF Report**

No

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been cleared or approved for diagnostic use by the U.S. Food and Drug Administration.

**CPT Code Information**

86008

---

**LOINC® Information**

| Test ID | Test Order Name           | Order LOINC® Value |
|---------|---------------------------|--------------------|
| FFWNC   | Walnut Component rJug r 3 | 81789-0            |

| Result ID | Test Result Name          | Result LOINC® Value |
|-----------|---------------------------|---------------------|
| Z5674     | Walnut Component rJug r 3 | 81789-0             |
| Z5675     | Class                     | 103118-6            |