



Test Definition: FPF12

Prothrombin Fragment 1+2 MoAb

Overview

Method Name

Enzyme-Linked Immunosorbent Assay

NY State Available

Yes

Specimen

Specimen Type

Plasma Na Cit

Specimen Required

Patient Preparation: Do not draw from an arm with a heparin lock or heparinized catheter.

Specimen Type: Citrated plasma

Collection Container/Tube: Light-blue top (citrate)

Specimen Volume: 2 mL

Collection Instructions:

1. Draw blood in a light blue-top (Sodium citrate) tube(s).
2. Centrifuge and aliquot 2 mL citrated plasma into a plastic vial.
3. Send frozen.

Note: To avoid delays in turnaround time when requesting multiple tests on frozen samples, **please submit separate frozen specimens for each test requested.**

Specimen Minimum Volume

1 mL (Note: This volume does not allow for repeat testing.)

Reject Due To

Specimens other than Plasma Anticoagulants other than sodium citrate	Reject
--	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	365 days	

Clinical & Interpretive**Reference Values**

<326 pmol/L

Pre-analytical conditions such as a difficult draw may spuriously increase test results.

Performance**PDF Report**

No

Day(s) Performed

Monday, Thursday

Report Available

5 to 13 days

Performing Laboratory Location

Esoterix Coagulation

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

83520

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPF12	Prothrombin Fragment 1+2 MoAb	27824-2

Result ID	Test Result Name	Result LOINC® Value
FPF12	Prothrombin Fragment 1+2 MoAb	27824-2