

Overview

Method Name

Quantitative Spectrophotometry

NY State Available

Yes

Specimen

Specimen Type

Body Fluid

Specimen Required

Specimen Type: Pericardial fluid

Collection Container/Tube: Leak-proof container

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions:

1. Centrifuge pericardial fluid specimen at ambient temperature.
2. Transfer 0.5 mL of fluid (supernatant) into a plastic vial and freeze.
3. Send frozen. **Specimen must remain frozen** until received at performing lab.
4. Indicate source on requisition.

Specimen Minimum Volume

0.2 mL

Reject Due To

Bronchoalveolar lavage (BAL)	Reject
Glass tubes	Reject
Turbid specimens	Reject
Whole blood	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Body Fluid	Frozen (preferred)	30 days	
	Refrigerated	7 days	

Clinical & Interpretive

Clinical Information

Refer to <https://ltd.aruplab.com/>

Reference Values

0-40 U/L

Performance

PDF Report

No

Day(s) Performed

Tuesday, Thursday, Sunday

Report Available

1 to 8 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

CPT Code Information

84311

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
---------	-----------------	--------------------

Test Definition: FADDP

Adenosine Deaminase, Pericardial Fluid

FADDP	Adenosine Deaminase Pericardial Fld	49760-2
Result ID	Test Result Name	Result LOINC® Value
Z5894	Adenosine Deaminase Pericardial Fld	49760-2