

Test Definition: FADDP

Adenosine Deaminase, Pericardial Fluid

Overview

Method Name

Quantitative Spectrophotometry

NY State Available

Yes

Specimen

Specimen Type

Body Fluid

Specimen Required

Specimen Type: Pericardial Fluid

Sources: Pericardial Fluid

Container/Tube: Standard Transport Tube

Specimen Volume: 0.5 mL

Collection Instructions: Collect Pericardial Fluid in a leak-proof container. Centrifuge specimen at room temperature,

transfer 0.5 mL pericardial fluid to plastic vial and freeze.

Note: Specimen must remain frozen until received at performing lab.

Specimen Minimum Volume

0.2 mL

Reject Due To

Hemolysis	NA NA	
Lipemia	NA	
Icterus	NA	
Other	Whole blood, Bronchoalveolar lavage (BAL) specimens, Turbid specimen	

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Body Fluid	Frozen (preferred)	30 days	
	Refrigerated	7 days	

Clinical & Interpretive



Test Definition: FADDP

Adenosine Deaminase, Pericardial Fluid

Reference Values

0-40 U/L

Performance

PDF Report

No

Day(s) Performed

Sunday, Tuesday, Thursday

Report Available

1 to 8 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

CPT Code Information

84311

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FADDP	Adenosine Deaminase Pericardial Fld	49760-2

R	Result ID	Test Result Name	Result LOINC® Value
Z	25894	Adenosine Deaminase Pericardial Fld	49760-2