



# Test Definition: FASPG

Aspergillus IgG Precipitins Panel

## Overview

### Method Name

Gel Diffusion

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

**Collection Container/Tube:**

**Preferred:** Red top

**Acceptable:** Serum gel

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL Serum

**Collection Instructions:**

1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerate.

### Specimen Minimum Volume

Serum: 0.5 mL

### Reject Due To

Gross hemolysis	Reject
Gross lipemia	Reject
Gross Icterus	Reject

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

## Clinical & Interpretive

### Clinical Information

Refer to [www.eurofins-viracor.com/test-menu/](http://www.eurofins-viracor.com/test-menu/)

### Reference Values

Negative

### Interpretation

The gel diffusion method was used to test this patients serum for the presence of precipitating antibodies (IgG) to the antigens indicated. These antibodies are serological markers for exposure and immunological sensitization. The clinical significance varies, depending on the history and symptoms.

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

7 to 11 days

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

### CPT Code Information

86331 x 6

---

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FASPG	Aspergillus IgG Precipitins Panel	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z5917	Aspergillus amstelodami/glaucus	Not Provided
Z5918	Aspergillus flavus	Not Provided
Z5919	Aspergillus fumigatus Mix	Not Provided
Z5920	Aspergillus nidulans	Not Provided
Z5921	Aspergillus niger	Not Provided
Z5922	Aspergillus versicolor	Not Provided