



Test Definition: FCDSU

Comprehensive Drug Screen, Umbilical Cord
Tissue

Overview

Profile Information

Test Id	Reporting Name	Available Separately	Always Performed
FCDU1	Opiates, UC	No	Yes
FCDU2	Benzodiazepines, UC	No	Yes
FCDU3	Acetyl Fentanyl, UC	No	Yes
FCDU4	Novel Psychoactives, UC	No	Yes
FCDU5	Amphetamines, UC	No	Yes
FCDU6	Cocaine and Metabolites, UC	No	Yes
FCDU7	Cannabinoids, UC	No	Yes
FCDU8	Dextro/Levo Methorphan, UC	No	Yes
FCDU9	Muscle Relaxants, UC	No	Yes
FCDUA	Phencyclidine, UC	No	Yes
FCDUB	Barbiturates, UC	No	Yes
FCDUE	Gabapentin, UC	No	Yes
FCDUD	Zolpidem, UC	No	Yes

Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
FCDV1	Acetyl Fentanyl Conf (Qual), UC	No	No
FCDV2	Amphetamines Conf (Qual), UC	No	No
FCDV3	Barbiturates Conf (Qual), UC	No	No
FCDV4	Benzodiazepines Conf (Qual), UC	No	No
FCDV5	Buprenorphine/Metab Conf (Qual), UC	No	No
FCDV6	Cannabinoids Conf (Qual), UC	No	No
FCDV7	Carisoprodol/Metab Conf (Qual), UC	No	No
FCDV8	Cocaine/Metab Conf (Qual), UC	No	No
FCDV9	Dextro/Levo Methorphan	No	No

	Conf Qual, UC		
FCDW1	Fentanyl/Metab Conf (Qual), UC	No	No
FCDW2	Meperidine/Metab Conf (Qual), UC	No	No
FCDW3	Methadone/Metab Conf (Qual), UC	No	No
FCDW4	Novel Psychoactive Conf (Qual), UC	No	No
FCDW5	Opiates-Free Unconj Conf (Qual), UC	No	No
FCDW6	Phencyclidine Conf (Qual), UC	No	No
FCDW7	Tapentadol Conf (Qual), UC	No	No
FCDW8	Tramadol/Metab Conf (Qual), UC	No	No
FCDW9	Zolpidem Conf (Qual), UC	No	No
FCDX1	Gabapentin Conf (Qual), UC	No	No
FCDX2	Mitragynine Conf (Qual), UC	No	No
FCDX3	Xylazine Conf (Qual), UC	No	No

Testing Algorithm

If appropriate, confirmation testing will be performed at no additional charge.

Method Name

High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

NY State Available

Yes

Specimen

Specimen Type

Tissue

Specimen Required

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Specimen Type: Umbilical cord tissue

Container/Tube: Plastic, preservative-free container

Specimen Volume: 10 g

Collection Instructions:

1. Collect at least 6 inches of umbilical cord (approximately the length of an adult hand).
2. Drain and discard any blood.
3. Rinse the exterior of the cord segment with normal saline or sterile water.
4. Pat the cord dry and place in a preservative-free plastic container for transport.
5. Ship refrigerated.

Specimen Minimum Volume

See Specimen Required

Reject Due To

Cords soaking in blood or other fluid	Reject
Formalin fixed	Reject
Tissue that is obviously decomposed	Reject
Samples received without cold pack	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Tissue	Refrigerated (preferred)	7 days	
	Frozen	14 days	

Clinical & Interpretive

Clinical Information

Refer to www.nmslabs.com/test-catalog

Reference Values

Reporting limit(s) determined each analysis.

None Detected

Performance**PDF Report**

No

Day(s) Performed

Monday through Sunday

Report Available

3 to 7 days

Performing Laboratory Location

NMS Labs

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80307

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCDSU	Comprehensive Drug Screen, UC	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6187	Zolpidem	97310-7
Z6185	Phenobarbital	97300-8
Z6186	Butalbital	97224-0
Z6184	Phencyclidine	97298-4
Z6182	Meprobamate	97270-3
Z6183	Carisoprodol	97226-5
Z6302	Xylazine	Not Provided

Z6181	Dextro/Levo Methorphan	97588-8
Z6179	Delta-9 Carboxy THC	97584-7
Z6180	Delta-9 THC	97586-2
Z6176	Benzoylcegonine	97222-4
Z6177	Cocaethylene	97232-3
Z6178	Cocaine	97234-9
Z6171	Amphetamine	97220-8
Z6172	MDA	97262-0
Z6173	MDEA	97264-6
Z6174	MDMA	97266-1
Z6175	Methamphetamine	97274-5
Z6168	alpha-PVP	97216-6
Z6169	Ethylone	97246-3
Z6170	Methylone	97276-0
Z6298	Mitragynine	Not Provided
Z6167	Acetyl Fentanyl	97214-1
Z6152	7-Amino Clonazepam	97212-5
Z6153	Alprazolam	97218-2
Z6154	Chlordiazepoxide	97228-1
Z6155	Clonazepam	97230-7
Z6156	Desalkylflurazepam	97238-0
Z6157	Diazepam	97240-6
Z6158	Flunitrazepam	97250-5
Z6159	Flurazepam	97252-1
Z6160	Hydroxytriazolam	97258-8
Z6161	Lorazepam	97260-4
Z6162	Midazolam	97278-6
Z6163	Nordiazepam	97282-8
Z6164	Oxazepam	97322-2
Z6165	Temazepam	97304-0
Z6166	Triazolam	97308-1
Z6131	6-Monoacetylmorphine-Free	97210-9
Z6132	Codeine-Free	97236-4
Z6133	Dihydrocodeine/Hydrocodol-Free	97242-2
Z6134	Hydrocodone-Free	97254-7
Z6135	Hydromorphone-Free	97256-2
Z6136	Morphine-Free	97280-2
Z6137	Norhydrocodone	97286-9
Z6138	Noroxycodone	97290-1
Z6139	Oxycodone-Free	97295-0
Z6140	Oxymorphone-Free	97296-8
Z6141	EDDP	97244-8

Z6142	Methadone	97272-9
Z6143	Fentanyl	97248-9
Z6144	Norfentanyl	97284-4
Z6145	Meperidine	97268-7
Z6146	Normeperidine	97288-5
Z6147	Norbuprenorphine-Free	97590-4
Z6148	Buprenorphine-Free	97582-1
Z6149	O-Desmethyltramadol	97292-7
Z6150	Tramadol	97306-5
Z6151	Tapentadol	97302-4
Z6299	Gabapentin	Not Provided