

## Overview

### Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
FA13R	Factor 13 1:1 Mix	No	No

### Testing Algorithm

If the factor 13 qualitative result states "Lysis", then factor 13 1:1 mix will be performed at an additional charge.

### Method Name

Qualitative Solubility Assay

### NY State Available

Yes

## Specimen

### Specimen Type

Plasma Na Cit

### Specimen Required

**Collection Container/Tube:** Light-blue top (3.2% sodium citrate)

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 2 mL Platelet-poor plasma

#### Collection Instructions:

1. Centrifuge immediately and aliquot 2 mL platelet-poor plasma in a plastic vial.
2. Send frozen.

**Additional Information:** Critical frozen. Separate specimens must be submitted when multiple tests are ordered.

### Specimen Minimum Volume

Platelet-poor plasma: 1 mL

### Reject Due To

Hemolysis	Reject
Clotted specimens	Reject
Serum	Reject
EDTA plasma	Reject

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

**Clinical & Interpretive****Clinical Information**

Refer to <https://ltd.aruplab.com/>.

**Reference Values**

Factor XIII, Qualitative: No Lysis

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

4 to 5 days

**Performing Laboratory Location**

ARUP Laboratories

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**CPT Code Information**

85291-Factor 13, Qualitative

85291-Factor 13, 1:1 Mix (if appropriate)

---

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FA13Q	Factor 13 Qual, Reflex to 1:1 Mix	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6312	Factor XIII, Qualitative	3241-7